

FLUSSO™

# By Pass Adapter

To maintain a closed ventilator circuit



Made in Canada

Flusso™ and Swing Valve Technology™ are trademarks used under licence by  
McArthur Medical Sales Inc.

US and International Patents Pending



## Flusso™ Highlights:

- Proprietary valve with Swing Valve Technology™
- Clear housing with colored valve permits visualization of function of Swing Valve Technology™
- Integrated swivel at patient connection to reduce torque on endotracheal tube
- Tethered port cap keeps ports clean and protects swing valve
- Minimal dead space (~6cc)
- Accommodates standard Heat and Moisture Exchangers (HME), resuscitation bag and breathing circuits
- Standard ISO connections

Flusso™ is a by pass adapter that can be used with passive or active humidity systems to prevent the disconnection of the patient from the mechanical ventilator during a circuit change or disconnection of patient for transport.



**FLUSSO™**  
BY PASS ADAPTER

[www.flussobypass.com](http://www.flussobypass.com)

To Order: Product Code 420-1978 25/case

*MMSI*  
McARTHUR  
MEDICAL SALES INC.

McArthur Medical Sales Inc  
1846 5th Concession Rd W  
Rockton, ON L0R 1X0

1-800-996-6674  
mmsi@mcarthurmedical.com  
www.mcarthurmedical.com



SMI1-420-1978