

medical scales

handrail scale

model 250-10-2



model 250-10-4



RICE LAKE
HEALTH SOLUTIONS

To be the best by every measure[®]

www.ricelake.com/health

handrail scale

THE DIGITAL HANDRAIL SCALE by Rice Lake is designed to weigh individuals who need extra support. The non-skid platform and side rails give extra stability and safety.

Advanced movement compensation technology adjusts for involuntary movement occurring when weighing unsteady individuals. Two heavy-duty wheels located at the rear of the base make maneuvering the handrail scale easy in any clinic or office.

The Digital Handrail Scale with hinged seat adds the unique ability for patients to be weighed in either a standing or sitting position.

	RL250-10-2	RL250-10-4
Capacity	1000lb (450kg)	1000lb (450kg)
Graduations	0.2lb (0.1kg)	0.2lb (0.1kg)



250-10 Series Standard Features



- Platform dimensions: 25 in W x 20 in L x 3 in H (64 cm W x 51 cm L x 8 cm H)
- Handrails for extra support and safety
- Integrated wheels for ease of portability
- Motiontrap™ (movement compensation technology)
- Large 1 inch LCD display (2.5 cm)
- Low-battery indicator/automatic power-off
- BMI (body mass index) function
- Powered by 6 "AA" batteries (included)
- Hold function
- Unit of measure: lb only, kg only, lb/kg
- EMR capable, USB
- Updatable firmware

RL250-10-2 Bariatric/Handrail

- Capacity: 1000lb x 0.2lb (450 kg x 0.1 kg)

RL250-10-4 Bariatric/Handrail with Chair Seat

- Capacity: 1000lb x 0.2lb (450 kg x 0.1 kg)
- Folding seat

Options/Accessories

- AC Adaptor

Warranty

Two-year limited warranty

Your Rice Lake Weighing Systems distributor is:



1846 5th Concession W., P.O. Box 7, Rockton, Ontario L0R 1X0
Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142
www.mcarthurmedical.com • mmsi@mcarthurmedical.com



To be the best by every measure®

230 W. Coleman St. • Rice Lake, WI 54868 • USA
TEL: 715-234-9171 • FAX: 715-234-6967 • www.ricelake.com

An ISO 9001 Registered Company © 2012 Rice Lake Weighing Systems PN 106426 11/12
Specifications subject to change without notice