



ET ADHESIVE TAPE™

Quick, effective,
securement in a
ONE-PIECE,
easy step



Three types of Marpac ET Tape are available:

- 301** – Latex-free silk
- 302** – Cloth featuring an extremely durable adhesive*
- 304** – Latex-free cloth featuring zinc oxide, is waterproof and hypoallergenic

“The Advanced Cardiac Life Support (ACLS) is recommending a manufactured product in all intubation trays.”

* CAUTION: The 302 product contains natural rubber latex which may cause allergic reactions.

New ET Tape from Marpac takes less time to apply than conventional cloth tape, and secures better, which reduces the chance of the ET tube becoming dislodged or extubation. The skin friendly 304 ET securement device includes zinc oxide, is waterproof and is hypoallergenic for the patient’s comfort.

ET Tape is individually packaged and has a unique one-piece design for quick application. It’s made with a soft cotton/foam



neck band and medical-grade adhesive tape. One end is enclosed in a rigid, disposable sleeve, which allows it to slide easily under the patient’s neck, and the scored backing is easy to peel off, even with protective gloves.



1846 5th Concession W., P.O. Box 7, Rockton, Ontario L0R 1X0
Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142
www.mcarthurmedical.com • mmsi@mcarthurmedical.com