

# McTree Tube Support™



The **McTree Tube Support™** works in conjunction with both adult and pediatric corrugated breathing tubes. It may be placed either above or below the surgical table and supports tubes in all positions. Assists in the reduction of tension on the endotracheal tube or face mask by helping to prevent kinks in the tube, accidental extubation and/or movement of the endotracheal tube. Made of surgical grade stainless steel, it is durable and can be easily cleaned.

Product Code: 420-1947



**Manages Adult  
and Pediatric  
Corrugated  
Breathing Tubes**

**Surgical Grade  
Stainless Steel**

**Latex Free**

**Reusable**

**Durable**

McArthur Medical  
Sales Inc.  
1846 5th Concession W.,  
Rockton, Ontario L0R 1X0

[www.mcarthurmedical.com](http://www.mcarthurmedical.com)

1-800-996-6674