

NeoTree™ Tube Support



The **NeoTree™ Tube Support** works in conjunction with neonatal corrugated breathing tubes. Reversible to secure tubes for term and pre-term neonates and supports tubes in all positions. Assists in the reduction of tension on the endotracheal tube or face

mask by helping to prevent kinks in the tube, accidental extubation and/or movement of the endotracheal tube. Made of lightweight aluminum with a baked-on powder coat finish, durable and can be easily cleaned.

Product Code: 420-1948



**Reversible to
Secure Term and
Pre-term Neonatal
Corrugated
Breathing Tubes**

**Lightweight
Aluminum**

**Easy to Clean
Powder Coat
Finish**

Reusable

Durable

McArthur Medical
Sales Inc.
1846 5th Concession W.,
Rockton, Ontario L0R 1X0

www.mcarthurmedical.com

1-800-996-6674