

## neonatal scale



model 650-10-1



**RICE LAKE®**  
HEALTH SOLUTIONS

To be the best by every measure®

[www.ricelake.com/health](http://www.ricelake.com/health)

## neonatal scale

### Accurate Weighing for Active Newborns and Infants

RICE LAKE'S RL-650 SERIES uses Motiontrap™ motion sensing technology to precisely weigh active babies. By compensating for involuntary movements, Rice Lake's neonatal scale is able to weigh newborns and infants up to 33 pounds. A two-sided removable tray keeps babies safe while being weighed and a built-in measuring tape ensures quick and easy measurements. Weight can be displayed in pounds/ounces or kilograms/grams. An optional stainless steel cart and AC adapter provide portability in any hospital or clinic setting.

<b>Capacity</b>	0-10lb (0-5 kg)	10-33lb (5-15 kg)
<b>Graduations</b>	0.05 oz (1 g)	0.1 oz (2 g)



Optional stainless steel cart

### 650-10-1 Standard Features



- Capacity: 0-10lb (0-5 kg) X 0.05 oz (1 g); 10-33lb (5-15 kg) X 0.1 oz (2 g)
- Tray dimensions: 11 x 23.5 x 3 in
- Two-sided removable tray with built-in measuring tape (up to 23.5 in x 1/8 in or 60 cm x 1/4 cm)
- Motiontrap (movement compensation technology)
- 0.75 in LCD display
- Reweigh function
- Low battery indication
- Auto shutoff
- Powered by (6) Li-FeS2 batteries or optional AC adapter

### Options/Accessories

- AC adapter
- Stainless steel cart

### Warranty

- Two-year limited warranty

### Approvals



1846 5th Concession W., P.O. Box 7, Rockton, Ontario L0R 1X0  
Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142  
www.mcarthurmedical.com • mmsi@mcArthurmedical.com



To be the best by every measure®

230 W. Coleman St. • Rice Lake, WI 54868 • USA  
TEL: 715-234-9171 • FAX: 715-234-6967 • [www.ricelake.com](http://www.ricelake.com)

An ISO 9001 registered company © 2015 Rice Lake Weighing Systems PN 106420 4/15  
Specifications subject to change without notice