

LEGAL
FOR TRADE

medical scales

portable
fitness
scale



model 140-10-7N



Measurement
Canada
Approved



RICE LAKE
HEALTH SOLUTIONS

To be the best by every measure®

www.ricelake.com/health

portable **fitness** scale



THE 140-10-7N PORTABLE FITNESS SCALE is an NTEP-Certified legal-for-trade device that provides the highest possible accuracy when weight is a critical factor. In wrestling programs, weight-loss centers and even health clubs, weight can determine membership fees or classification levels. For that reason a legal-for-trade device must be used.

This portable legal-for-trade device offers a full 600 lb capacity, making it the perfect tool for a broad range of weighing needs. Mount your display on a wall or place on a table for fast and convenient use. Weighing only 20 lbs, this scale base and indicator package can easily be carried from event to event. Weight is displayed in both standard and metric values, accurate to within 0.1lb. An onboard Body Mass Index (BMI) calculator helps save time and increase health awareness. Operate up to 75 hours with internal rechargeable battery or plug in at any location. An optional carrying case makes for easy transportation of this lightweight and compact scale.

140-10-7N Standard Features

- Platform dimensions: 14 in x 14.5 in x 3 in
- Large 1 inch Liquid Crystal Display (LCD)
- AC adaptor and internal rechargeable battery
- RS-232 output
- 6.5 ft cable between scale base and indicator
- Switch between lb and kg

Options/Accessories

- Transport/carrying case

Warranty

- Two-year limited warranty

Approvals



Measurement
Canada
Approved



Capacity	600 lb (270 kg)
Graduations	0.1 lb (0.05 kg)

weight loss



fitness



athletics



Your Rice Lake Weighing Systems distributor is:



Canadian Distributor:
 McArthur Medical Sales Inc.
 1846 5th Concession West, Rockton, ON L0R 1X0
 T+ 800.996.6674 519.622.4030 F+ 519.622.1142
 www.mcarthurmedical.com mmsi@mcarthurmedical.com

RICE LAKE
HEALTH SOLUTIONS

To be the best by every measure®

230 W. Coleman St. • Rice Lake, WI 54868 • USA
 TEL: 715-234-9171 • FAX: 715-234-6967 • www.ricelake.com

An ISO 9001 Registered Company © 2010 Rice Lake Weighing Systems PN 114089 4/10
 Specifications subject to change without notice