



# Precision Syringes



## For peace of mind regarding the accuracy of your spirometer

Vitalograph precision syringes provide quick and easy accuracy checks with reliability and fast verification.

Our precision syringes are available as a simple 3-L or a valved 1-L type. The 3-L precision syringe is for all spirometers as required by ISO 26782:2009. It can be used with flow or volume sensor systems and can be used with any manufacturers' equipment. A range of adaptors are available. The 1-L Precision Syringe is for linearity checking volume displacement spirometers.

- Certified accuracy better than 0.5%
- End-to-end volume, no lining up marks or grooves
- Traceable to international standards
- Unique serial number and history file for each syringe



1846 5th Concession W., P.O. Box 7, Rockton, Ontario L0R 1X0  
Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142  
[www.mcarthurmedical.com](http://www.mcarthurmedical.com) • [mmsi@carthurmedical.com](mailto:mmsi@carthurmedical.com)