



the difference is clear

Vitalograph copd-6

Revolutionary COPD screener

- ❖ Identifies those at risk of COPD at the pre-symptomatic stage to allow early medical intervention and facilitate better clinical outcomes
- ❖ Screens out those who's FEV₁ is normal, who therefore do not have COPD, without the risk of false negatives
- ❖ Diagnostic spirometry resources can be focused on those most likely to be diagnosed with COPD
- ❖ Displays the key parameters for clinical interpretation including obstructive index and COPD classification (stage I – IV)
- ❖ Displays FEV₁, FEV₆, ratio and % predicted, obstructive index, COPD classification and lung age
- ❖ Detachable flowhead makes cleaning simple
- ❖ Easy to use by non-respiratory specialists
- ❖ Built-in quality of blow indicator
- ❖ Designed for multi-subject use
- ❖ Large, easy to read display

Vitalograph is the world's leading provider of outstanding respiratory diagnostic, screening and monitoring devices. With a heritage of excellence spanning half a century, Vitalograph continue to make pioneering contributions to effective respiratory care and enhanced quality of life.





the difference is clear



Technical Specification:

Product:
Vitalograph copd-6

Model Number:
4000

Parameters Displayed:
FEV₁/FEV₆/ratio and their predicted values

Quality of Blow Indicator:
Yes

Obstructive Index:
Mild; Moderate; Severe (FEV₁% Predicted)

COPD Classification:
Normal (not COPD); Stage I, II, III & IV (FEV₁%Predicted)

Lung Age Indicator:
Yes

Accuracy:
Better than +/-3%

Range:
0 - 9.99 L BTPS

Sensor:
Stator/rotor

Flow Impedance:
Better than 0.15 kPa/L/s at 14 L/s

Power Supply:
AAA Batteries

Display:
Custom Liquid Crystal Display

Size:
113 x 63 x 48 mm

Weight:
55g

Performance Standards:
ISO EN 13826:2003

ATS/ERS Guidelines:
2005

Operating Temperature:
17 -37 °C

Safety Standards:
EN ISO 60601-1:2005 -1 -2 -4

Medical Safety Standard:
Medical Devices Directive
93/42/EEC1993

Designed & Manufactured Under:
ISO 13485:2003, FDA 21CFR820,
CMDR & JPAL

Ordering Info:

Catalogue No.

- | | |
|-------|---------------------------------|
| 40200 | Vitalograph copd-6 |
| 20242 | SafeTway® mouthpieces (200) |
| 20980 | Mini SafeTway® mouthpieces (50) |
| 20303 | Noseclips (200) |

USB VERSION AVAILABLE



1846 5th Concession W., P.O. Box 7, Rockton, Ontario L0R 1X0
Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142
www.mcarthurmedical.com • mmsi@mcarthurmedical.com

www.vitalograph.co.uk

