


medical scales

RL-MPS

Healthweigh™
mechanical physician scale



RICE LAKE™
WEIGHING SYSTEMS

Commitment Beyond Measurement®

www.ricelake.com

RL-MPS Healthweigh mechanical physician scale

THE RICE LAKE MECHANICAL Physician Scale features a dual die-cast beam for weighing up to 440 lb (200 kg). The weight can be read from either the front or back. The scale base is designed with a removable slip-resistant plastic base cover for easy cleaning.

The standard height rod measures 24 (60 cm) to 84 inches (212 cm). Durable, accurate and easy to operate, move and store, the Rice Lake Healthweigh Mechanical Physician Scale is ideal for use in clinics and hospitals.

Capacity	440 lb (200 kg)
Graduations	0.25 lb (0.1 kg)


Standard Features

- Dual read capacity 440 lb x 0.25 lb (200 kg x 0.1 kg)
- 10.75" x 14.75" platform
- Retractable aluminum height rod
- Sturdy enameled steel body
- Die-cast beam
- Rugged scale understructure
- Removable slip-resistant plastic case cover
- Easy operation with two sliding weights
- Rear wheels for ease of portability

Shipping information

Shipping Weight	46 lb
Shipping Dimension	60" x 21" x 7"

Heavy-duty wheels allow easy portability of the scale.



Multifunctional, height and weight are indicated in a one task operation. Weigh information can be read from the front and back.



Master Distributor - Canada:

McArthur Medical Sales Inc.
 1846 5th Concession West
 Rockton, ON L0R 1X0
 T+ 800.996.6674 F+ 519.622.1142
mmsi@mcarthurmedical.com