

digital **baby** scale



RICE LAKE
WEIGHING SYSTEMS

Commitment Beyond Measurement®

www.ricelake.com/health

digital **baby** scale

EVEN THE MOST ACTIVE BABIES can be weighed and measured accurately, safely and quickly with the RL-DBS. The baby can be placed on the scale, quickly removed, and then have its weight recalled for charting by pressing the *Recall* key. Rice Lake's battery or AC-powered Digital Baby Scale can be placed anywhere in the facility. The RL-DBS is sleek, extremely durable, and easy to clean and maintain.

Capacity	44 lb (20 kg)
Graduations	0.5 oz (10 g)

RL-DBS Standard Features

- Platform Dimensions: 11 in W x 23.5 in L
- Convenient lb/kg toggle key
- Large, easy-to-read 1 inch Liquid Crystal Display (LCD)
- Built-in measuring tape 0-22.5 inches (0-57 cm)
- *Recall* key to view last weight measurement
- Power: Four AA batteries (not included) or AC adaptor (included)
- Automatic shut-off ensures long battery life

Warranty

- Two-year limited warranty



Your Rice Lake Weighing Systems distributor is:



McArthur Medical Sales Inc.
1846 5th Concession West
Rockton, ON L0R 1X0
T+ 800.996.6674 519.622.4030
F+ 519.622.1142
www.mcarthurmedical.com
mmsi@mcarthurmedical.com



Commitment Beyond Measurement®

230 W. Coleman St. • Rice Lake, WI 54868 • USA
 TEL: 715-234-9171 • FAX: 715-234-6967 • www.ricelake.com

© 2009 Rice Lake Weighing Systems PN 110619 4/09
 Specifications subject to change without notice