

CALYPSO IS THE TEST WINNER

NOISE TEST SHOWS THAT CALYPSO IS THE QUIETEST BREASTPUMP*



* Visit www.ardomedical.com/silence-test
for the test conditions and detailed results

SWISS QUALITY 

 **ARDO**
Hospital Recommended

CALYPSO PROVES TO BE THE QUIETEST

MORE PEACE AND QUIET FOR MOTHER, FATHER AND BABY

Would you like more peace and quiet for you and your baby?

Ardo has the answer: the Calypso – either the double or single pump version – is the quietest of all¹ electric breastpumps. In comparison with many other commercially available pumps, it is only perceived to be half as loud², but it is still powerful.


Much-needed peace and quiet

The Calypso is quieter than all the other tested pumps. This is confirmed by the test report by ISO-certified, independent engineering firm Müller-BBM from Planegg near Munich, which refers to a validated noise level test of a range of electric breastpumps conducted in 2014. According to market research³, the level of noise made by a breastpump is an important feature.

In May 2014, Ardo medical AG commissioned the German measurement institute to test the sound power level of the most popular electric breastpumps. Along with the Calypso and the Calypso Double Plus, the institute investigated comparable breastpumps from other providers, too.

The result speaks for itself

The investigation by measurement institute Müller-BBM compares the noise levels of the top commercially available electric breastpumps. The measurements, which were taken in expression mode at the maximum vacuum setting in each case, show that the Calypso is quieter than all the other pumps tested: the Calypso single pump reaches a noise level of only 46.0 decibels, while the Calypso double pump reaches a mere 46.4 decibels – about the same as a radio set to a low volume. While in use, the Calypso is noticeably quieter than comparable products – and now the test results confirm this.

 Hear for yourself with our audio files:
www.ardomedical.com/silence-test

¹ Refers to all leading comparable breastpumps known to the company at the time of the test.

² The human ear cannot hear all frequencies. Specialists therefore assume that +10 dB is double the volume as heard by the human ear. Psychoacoustic effects are also taken into consideration in this calculation.

³ Breast Pumps Market, Copyright © 2013 Transparency Market Research


Quiet and powerful, yet pleasantly gentle

If a new mother wants to express milk regularly and for long periods of time without any restrictions, the Calypso breastpump is the most convincing choice: the Calypso is the quietest electric breastpump we know of. It therefore offers more peace and quiet for mother, father and baby.


A comparison of the breastpumps tested

You can find information on the validated test conditions and detailed test results at www.ardomedical.com/silence-test

Single pumps



Double pumps



Our quiet breastpumps:

CALYPSO

ELECTRIC BREASTPUMP



CALYPSO DOUBLE PLUS

DOUBLE ELECTRIC BREASTPUMP



*Happy Mothers for
Happy Babies*



McARTHUR MEDICAL SALES INC.

TF / 800.996.6674 P / 519.622.4030

mmsi@mcarthurmedical.com | www.mcarthurmedical.com

1846 Concession 5 W | P.O. Box 7 | Rockton, ON | L0R 1X0

 **ARDO**
Hospital Recommended