



AccuPAP™

Lung Expansion Therapy in the Palm of Your Hand

- Provides positive airway pressure for the prevention and treatment of atelectasis.
- Optimizes delivery of aerosolized medications.
- PEEP independent of flow.
- Four PEEP Settings: 5/10/15/20cmH₂O.

AccuPAP™ - Positive Airway Pressure Device,

| Part No. | Product Description | Case/Qty. |
|----------|--|-----------|
| 313-6000 | AccuPAP with Mouthpiece and 7' Tubing | 10 |
| 313-6001 | AccuPAP with Mouthpiece, 7' Tubing, and DPM™ (Disposable Pressure Manometer) | 10 |



DPM™ - Disposable Pressure Manometer

| Part No. | Product Description | Case/Qty. |
|----------|-------------------------------------|-----------|
| 313-5558 | DPM (Disposable Pressure Manometer) | 10 |

