




BiTrac™ Full Face & Nasal Masks

Fixed and Select Options to Suit a Variety of Patient Needs

- The BiTrac Full Face and Nasal Masks come in a variety of sizes.
- Both product lines are available with fixed elbows in place.
- Our BiTrac Select™ Full Face Masks come with interchangeable elbows and a variety of cushion styles, to suit multiple patient needs.



BiTrac™ Nasal with fixed Anti-Asphyxia Elbow

Part No.	Sizes	Case/Qty.
313-9013W	Adult Small	10
313-9014W	Adult Medium	10
313-9015W	Adult Large	10
313-9016W	Adult Extra Large	10

BiTrac™ Full Face with Fixed Anti-Asphyxia Elbow

Sizes	Case/Qty.
313-9001W	Adult Small
313-9002W	Adult Medium
313-9003W	Adult Large
313-9004W	Adult Extra Large

BiTrac™ Full Face with Fixed Standard Elbow

Part No.	Case/Qty.
313-9009W	Adult Small
313-9010W	Adult Medium
313-9011W	Adult Large
313-9012W	Adult Extra Large

BiTrac Select™ Full Face with Interchangeable Anti-Asphyxia with leak Elbow

Part No.	Sizes	Case/Qty.
313-9511	Adult Small	10
313-9512	Adult Medium	10
313-9513	Adult Large	10
313-9514	Adult Extra Large	10

BiTrac Select™ Full Face with Interchangeable Standard Elbow

Part No.	Sizes	Case/Qty.
313-9506	Adult Small	10
313-9507	Adult Medium	10
313-9508	Adult Large	10
313-9509	Adult Extra Large	10

BiTrac Select™ Full Face with Interchangeable Anti-Asphyxia Elbow

Part No.	Sizes	Case/Qty.
313-9501	Adult Small	10
313-9502	Adult Medium	10
313-9503	Adult Large	10
313-9504	Adult Extra Large	10

BiTrac Select™ Interchangeable Elbows

Part No.	Sizes	Case/Qty.
313-9593	ICE Anti-Asphyxia w/ Leak Elbow	10
313-9594	ICE Anti-Asphyxia Nebulizer Elbow	10
313-9595	ICE Standard Nebulizer Elbow	10
313-9596	ICE Anti-Asphyxia Bronchoscopy Elbow	10
313-9597	ICE Standard Bronchoscopy Elbow	10
313-9598	ICE Anti-Asphyxia Elbow	10
313-9599	ICE Standard Elbow	10

