

Pulmodyne®

GO-PAP™

Independent PEEP with a Low Flow Connection.

The GO-PAP is a disposable emergency CPAP device, that runs off of a low flow connection. It operates at 10LPM, a setting that all regulators can agree with! Like our O2-MAX™ this device delivers approximately 30% FiO₂, and provides the option of 5, 7.5, and 10cm of PEEP.

The PEEP settings are **independent of flow**, so you never have to worry about adjusting your flow to maintain your PEEP levels. This eliminates the need to read from confusing charts, or monitor manometers.

Just set it and forget it, and concentrate on caring for your patient.



The GO-PAP System includes:

GO-PAP Generator
Oxygen Supply Line
BiTrac ED™ Mask and Headstrap
OmniClip™

* Purchase with or without
Nebulizing Cannister

Part No.	Product Description	Case/Qty.
Go-PAP with included nebulization cannister		
313-4602N	GO-PAP with BiTrac ED™ mask size Adult Medium	10
313-4603N	GO-PAP with BiTrac ED™ mask size Adult Large	10
Go-PAP without nebulization cannister		
313-4602	GO-PAP with BiTrac ED™ mask size Adult Medium	10
313-4603	GO-PAP with BiTrac ED™ mask size Adult Large	10