



**NOW YOU CAN
ATTACH ALL YOUR
RESPIRATORY
DEVICES!**



TORTLE NICU MIDLINER

Precision Head Positioning System for micro-preemies & preemies.

Midliner Benefits:

**Supports Supine Midline Position/
Neutral head positioning**

- The AAP recommends supine midline position for micro-preemies during the first 72+ hours of life to help prevent intra-ventricular hemorrhage (IVH)

**Velcro front opening provides flexible
size adjustment**

- Simple and easy to apply which minimizes handling and stress on baby
- Provides easy access for scalp IVs and head ultrasounds
- Ensures a proper and comfortable fit without slippage
- Reduces intervention time for patients & staff

New design + materials

- Securely holds all nasal CPAP interfaces and nasal cannulas including RAM, as well as feeding tubes
- Includes additional straps to hold most bilishades
- X-ray friendly and compatible with ventilator

**Simplifies transition to all positions, including
supine midline, side lying and prone (log rolling)**

- No need to detach and re-apply respiratory or feeding apparatus
- Helps prevent and treat Dolichocephaly/Scaphocephaly & Head preference

Available in three sizes

SIZE	CASE#	HEAD CIRCUMFERENCE
ELBW	65911	17-22 cm
MICRO	65901	22-30 cm
PREEMIE	65903	30-36 cm

- 10 units per case
- Includes wash bag & cannula attachments
- Wash in residential washers & dryers on cool cycle

Supports AAP Back-to-Sleep Campaign

Respiratory Kits

TYPE	CASE#	COMPATIBLE
MIDLINE	65919	F&P, DRAGER, MINIFLOW, CARE FUSION
GOAL POST	65921	HUDSON, BIBI



The FDA-cleared Tortle NICU Midliner is a simple and efficient head repositioning system that assists staff in providing optimal care during IVH protocol & also to help prevent dolichocephaly & head preference.

Three Primary Positions

1- Supine Midline

First open front and pull straps up over the side of the head & gently close straps.

Cannula Keys supplied can be used to secure cannulas & feeding tubes.



Extra Velcro is provided.

2- Side Lying



Velcro supplied can be used to support rolls.

3- Prone



CE FDA-cleared



1846 5th Concession W., P.O. Box 7, Rockton, Ontario L0R 1X0
Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142
www.mcarthurmedical.com • mmsi@mcarthurmedical.com