


# Neonatal Transilluminator

Specially Designed for Special Patients



Finally, a neonatal transilluminator designed specifically with you and the neonate in mind. We've taken all the great features of the **Venoscope® II** and put them into a transilluminator that better serves the special needs of your neonatal patients. It still features the same powerful LED lights, but is now smaller and more portable, making it even easier to use with neonates. The cool lights are strong enough to illuminate the hardest to find veins without the risk of burning the baby. Creating happier babies makes happier nurses.

The **Neonatal Transilluminator** weighs just 6 ounces including the 3 "AA" batteries. It features a low battery indicator light that alerts the nurse to change the batteries.

## Features:

- **Easy to use**
- **Bright / Cool lights that won't burn the baby**
- **Compact / Portable**  
**Battery Operated**

## Venoscope Disposables

### Disposable Protective Covers

The Disposable Protective Covers are made of clear plastic, no latex, and are designed to cover the Venoscope® II so that blood or body fluids do not contaminate the light. The covers are packaged in a convenient tear away pack of 50 covers per pack. Simply tear one of the covers off the pack and slip the Venoscope® II into it. The cover is designed to fold over to cover the entire light. Pull the cover over the arms so that there is a clear field of vision between the arms when in use.



### Velstretch Straps

The Velstretch Straps allow the nurse to free up both hands for the venipuncture and to help stabilize the vein to prevent vein rolling.

Only Available from

### McArthur Medical Sales Inc.

1846 5<sup>th</sup> Concession West

Rockton, ON L0R 1X0

T+ 800.996.6674

F+ 519.622.1142

E+ [mmsi@mcarthurmedical.com](mailto:mmsi@mcarthurmedical.com)

W+ [www.mcarthurmedical.com](http://www.mcarthurmedical.com)