

Opticare® and Opticare Arthro

Eye drop dispensers



For many people
Opticare® and
Opticare Arthro
solve the difficulties
associated with
administering
eye drops

NEW

An award winning design, Opticare® and Opticare Arthro offer patients significant benefits

Patient independence

Opticare® enables the accurate self administration of eye drops by many patients who would otherwise find it impossible. This frees up valuable nursing time.

Comfortable and hygienic

The dispensers fit securely around the eye and the in situ bottle positioning dispels the fear of the dropper tip touching the eye. Touch contamination is prevented.

Accurate targeting of medication

The anatomically designed eye piece raises the upper lid when placed against the eye. There is also finger space in the eye piece for pulling down the lower lid.

Easy to use

Medication bottles can be loaded easily and a double squeeze mechanism provides a controlled delivery of eye drops using either fingers or a grip.

Opticare®

Better by design

McArthur Medical Sales Inc.
1846 5th Concession West, Rockton, Ontario L0R 1X0
Tel: 1-800-996-6674 (toll free) www.mcarthurmedical.com



The Opticare® Range

Take control of eyedrops

Opticare®

The universal model. Accommodates most 2.5ml to 15ml plastic dropper bottles. Cylindrical with a grip-sized body to increase user confidence.




Opticare Arthro 5

Extending arms leading to a rotating eyepiece. Particularly helpful for patients with arthritis and those with limited hand/shoulder action. Accommodates most 2.5ml and 5ml plastic dropper bottles.



Opticare Arthro 10

Same design specifications as the Arthro 5. Accommodates most 10ml and 15ml plastic dropper bottles.



*For more information or to place an order please contact
McArthur Medical at the address and contact number below*

Opticare®

Better by design

McArthur Medical Sales Inc.
1846 5th Concession West, Rockton, Ontario L0R 1X0
Tel: 1-800-996-6674 (toll free) www.mcarthurmedical.com

