

DRI Sleeper® Bedwetting Alarms



Drops of Urine
Urosensor™

Signal to Alarm Unit
(By wire / radio signal)

DRI Sleeper®
Alarm Unit

Alarm
sounds

Child wakes up &
trains brain to wake

DRI Sleeper® Eclipse



The unique DRI Sleeper® Eclipse features:

No wires between enuresis alarm and sensor

Long-term battery in the Urosensor™ of up to 3 years (not one month as in some other alarms)

Patented Urosensor™ technology. Unlike almost all other sensors it has no metal parts that can be corroded by urine, or irritate the skin.

The Urosensor™ is a moulded and completely sealed unit. Leakage and the potential for damage to the electronics is a problem with most other product.

However, with the Urosensor™ urine cannot damage the electronics. Because it is moulded and sealed it is quick and simple to clean for instant reuse.

The special electronics make the wireless Urosensor™ the thinnest bedwetting sensor of all.

High volume with a tamper-proof volume adjustment.

The Urosensor™ can transmit to more than one DRI Sleeper® alarm unit, and so there can be another receiver in a parent's or carer's room. Large detecting area to detect the urine.

Unit size: 80 x 40 x 20mm (3 x 1.5 x 0.75 inches)

DRI Sleeper® Excel



The DRI Sleeper® Excel features:

Max alerting sound (we use the sound that wakes the best).

Volume control (we set it to max - the control is internal so your child can't change it).

Child friendly colours.


Unit size: 60 x 40 x 15mm (2.25 x 1.5 x 0.5 inches)

DRI Sleeper® has the only independently audited user survey: 90% success rate


Safety electronics (see reverse)

DRI *Sleeper*[®] Safety Electronics

DRI *Sleeper*[®] Excel - Wireless Bedwetting Alarm



DRI *Sleeper*[®] Eclipse - Wireless Bedwetting Alarm



MMSI
McARTHUR
MEDICAL SALES INC.

DISTRIBUTED BY:

McARTHUR MEDICAL SALES INC.

TF / 800.996.6674 P / 519.622.4030

mmsi@mcarthurmedical.com | www.mcarthurmedical.com

1846 Concession 5 W | P.O. Box 7 | Rockton, ON | LOR 1X0

Developed by
Dr. Anthony Page
by Anzacare certified
to ISO 13485:2003