

Extended Dwell Catheter Product Offerings

Providing optimal outcomes for long-term
and intermediate therapies.



Table of Contents

Argon PICC/Midline Insertion Trays..... 2

Argon Safety MST Accessory Kits for PICC and Midline Placement..... 3

Argon Silicone Catheters..... 4

Argon L-Cath™ Polyurethane Catheters..... 5

Argon Neonatal/Pediatric PICC and Midline Catheters including BD™ Introducers..... 6

Argon Clinical Education..... 8

Argon PICC/Midline Insertion Trays

Argon PICC/Midline Insertion Trays are designed with high quality components to meet your clinical needs during the insertion process to improve your patient outcomes. The Neonatal Maximal Barrier Insertion tray meets the INS Standards, CDC Guidelines, National Patient Safety Goals, NANN Guidelines and the Joint Commission - Preventing Central Line - Associated Bloodstream Infections. The Neonatal Maximal Sterile Barrier PICC/Midline Insertion Tray is not made with natural rubber latex, DEHP or lint.

PICC and Midline Insertion Trays (Does not include catheter) (Box of 5)

CATALOG # DESCRIPTION

Neonatal/Infant PICC or Midline Maximal Barrier Insertion Tray

Contents:

384835	1 Non-Strike-Thru CSR Wrap	1 Forceps with Tip Protectors	1 Microbore Extension Set (DEHP free/lipid resistant/not made with natural rubber latex)
	1 Facemask with Eye Shield	4 4" x 4" Absorbent Pads	
	1 Gown	4 2" x 2" Absorbent Pads	
	1 Bouffant Cap	1 CHG Swabsticks (3 per pack)	2 Pre-Filled 10 ml syringes with preservative-free 0.9% sodium Chloride Solution
	1 Pair Sterile Gloves	1 Alcohol Swabsticks (3 per Pack)	
	2 Poly-lined Drapes	1 Iodophor PVP Swabsticks(3 per pack)	
	1 Fenestrated Drape	1 Skin Protectant Swabstick	
	1 Tape Measure	1 Sterile Normal Saline Wipe	
	1 Alcohol Pad	1 1 Roll Transpore Tape	
	1 Tourniquet	1 Deadender Cap	
	1 Scissors	2 Semipermeable Transparent Dressings	

Maximal Barrier Insertion Tray - Adult/Pediatric PICC or Midline

Contents:

384837	1 Non-Strike-Thru CSR Wrap	1 Forceps with Tip Protectors	1 Microbore Extension Set (DEHP free/lipid resistant/not made with natural rubber latex)
	1 Facemask with Eye Shield	4 4" x 4" Absorbent Pads	
	1 Isolation Mask	4 2" x 2" Absorbent Pads	
	1 Gown	2 3ml Chloraprep	2 Pre-Filled 10 ml syringes with preservative-free 0.9% sodium Chloride Solution
	1 Bouffant Cap	1 1ml Syringe with Safety Needle	
	1 Pair Sterile Gloves	1 Preprinted Lidocaine Label	
	1 Non-Strike-Thru Maximal Barrier Drape	1 Skin Protectant Swabstick	1 Statlock Catheter Securement Device
	1 Poly-Lined Drape	1 Sterile Normal Saline Wipe	
	1 Fenestrated Drape	1 Roll Transpore Tape	
	1 Tape Measure	4 Foam Strips	
	1 Tourniquet	1 Semipermeable Transparent Dressing	
	1 Scissors		

Argon Silicone Catheters

Providing comfort for your patient

Argon First PICC™ and Argon First Midcath™ Catheters are high-quality silicone catheters designed for safe and easy placement while providing comfort to the patient.

- **Universal application-breadth** of line offers access options for neonate to geriatric patient population
- **Soft silicone integrated extension set** enhances patient comfort and to provide increased mobility
- **Pre attached T-port extension** allows flushing of the catheter before, during and after insertion
- **Protective sleeve** minimizes accidental contamination of the catheter
- **Centimeter markings** allow easy measurement and catheter trimming prior to insertion
- **Trim tool** is designed to ensure clean and accurate trimming and provides a 90° cut consistent with INS standards and FDA recommendation
- Not made with natural rubber latex and DEHP free



PICCs (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384232	26 G	1.9F	50 cm	S/L First PICC	
384242	26 G	1.9F	50 cm	S/L First PICC	Bard StatLock®
384221	26 G	1.9F	50 cm	S/L First PICC	BD Introsyte-N™ Introducer
384101	22 G	2.8F	50 cm	S/L First PICC	
384141	22 G	2.8F	50 cm	S/L First PICC	Bard StatLock
384156	20 G	3F	65 cm	S/L First PICC	
384115	20 G	3F	65 cm	S/L First PICC	BD Introsyte™ AutoGuard™ Introducer
384157	18 G	4F	65 cm	S/L First PICC	
384158	16 G	5F	65 cm	S/L First PICC	
384171	18 G	4F	65 cm	S/L First PICC	BD Introsyte AutoGuard Introducer
384177	16 G	5F	65 cm	S/L First PICC	BD Introsyte AutoGuard Introducer
384159	16 G	5F	65 cm	D/L First PICC	
384183	16 G	5F	65 cm	D/L First PICC	BD Introsyte AutoGuard Introducer

PICCs (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384357	20 G	3F	20 cm	S/L First Midcath Midline	
384304	20 G	3F	20 cm	S/L First Midcath Midline	BD Introsyte AugoGuard Introducer
384358	18G	4F	20 cm	S/L First Midcath Midline	
384308	18 G	4F	20 cm	S/L First Midcath Midline	BD Introsyte AutoGuard Introducer

Argon L-Cath™ Polyurethane Catheters

Breadth of line and familiarity you were trained on


Polyurethane material provides strength and increased flow rates. Argon L-Cath™ PICC and Midline Catheters are designed for safe and easy placement.

- **Universal application-breadth** of line offers access options for neonate to geriatric patient population
- **Pre attached Y-site** connector allows flushing of catheter before, during, and after the insertion
- **Protective sleeve** minimizes accidental contamination of the catheter
- **Trim tool** is designed to ensure clean and accurate trimming and provides a 90° cut consistent with INS Standards and FDA recommendation
- **Centimeter markings** allow easy measurement and trimming prior to insertion
- Not made with natural rubber latex and DEHP Free

PICCs (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384516	28 G	1.2F	25 cm	S/L L-Cath™ PICC	U-Wing Splittable Introducer Needle
384539	26 G	1.9F	30 cm	S/L L-Cath™ PICC	
384466	20 G	2.6F	60 cm	D/L L-Cath™ PICC	
384566	20 G	2.6F	60 cm	D/L L-Cath™ PICC	BD Introsyte™ AutoGuard Introducer
384463	20 G	3F	60 cm	S/L L-Cath™ PICC	
384481	20 G	3F	60 cm	S/L L-Cath™ PICC	BD Introsyte™ AutoGuard Introducer
384464	18 G	3.5F	60 cm	S/L L-Cath™ PICC	
384467	18 G	3.5F	60 cm	D/L L-Cath™ PICC	
384490	18 G	3.5F	60 cm	S/L L-Cath™ PICC	BD Introsyte™ AutoGuard Introducer
384575	18 G	3.5F	60 cm	D/L L-Cath™ PICC	BD Introsyte™ AutoGuard Introducer
384465	16 G	5F	60 cm	S/L L-Cath™ PICC	
384468	16 G	5F	60 cm	D/L L-Cath™ PICC	
384499	16 G	5F	60 cm	S/L L-Cath™ PICC	BD Introsyte™ AutoGuard Introducer
384584	16 G	5F	60 cm	D/L L-Cath™ PICC	BD Introsyte™ AutoGuard Introducer

Midline Catheters (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384620	20 G	3F	20 cm	S/L L-Cath Midline	
384631	20 G	3F	20 cm	S/L L-Cath Midline	BD Introsyte™ AutoGuard Introducer
384621	18 G	3.5F	20 cm	S/L L-Cath Midline	
384635	18 G	3.5F	20 cm	S/L L-Cath Midline	BD Introsyte™ AutoGuard Introducer

Complete Kits (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384518	28 G	1.2F	25 cm	L-Cath PICC Complete Kit	U-Wing Splittable Needle Introducer and 384835



Argon Neonatal/Pediatric PICC and Midline Catheters

Argon Neonatal/Pediatric PICC and Midline catheters are designed specifically for this patient population. Catheters come in a choice of silicone material or polyurethane material. The option for safety and non-safety introducers is designed to enhance first-stick proficiency. The BD™ Introducer Instaflash™ Needle Technology provides rapid visualization of blood return along the catheter shaft and is designed to enhance first-stick success in highly vulnerable patients.

Silicone and Polyurethane

Vascular access products you can trust

PICC and Midline Catheters including BD™ Introducers with Instaflash™ Needle technology

- **Soft silicone integrated extension set** designed to enhance patient comfort and to provide increased mobility
- **Pre-attached Y-site** on polyurethane catheters allows for flushing of catheter before, during, and after the insertion
- **Protective sleeve** minimizes accidental contamination of the catheter during placement
- **Trim tool** ensures clean and accurate trimming and provide a 90° cut consistent with INS Standards and FDA recommendations
- **Centimeter markings** allows easy measurement and catheter trimming prior to insertion
- Not made with natural rubber latex and DEHP free
- **Pre-attached T-Port extension** on silicone catheters allows flushing of the catheter before, during, and after insertion

Silicone and Polyurethane

Introducers (Box of 10)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384061	28 G	1.2F		U-Wing Splittable Introducer Needle	
384062	26 G	1.9F		FirstStic-N Introducer	

Silicone PICCs (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384232	26 G	1.9F	50 cm	S/L First PICC™	
384242	26 G	1.9F	50 cm	S/L First PICC™	Bard StatLock®
384221	26 G	1.9F	50 cm	S/L First PICC™	BD Introsyte-N™ Introducer
384222	26 G	1.9F	50 cm	S/L First PICC™	Bard StatLock® and Argon FirstStic N-Introducer
384101	22 G	2.8F	50 cm	S/L First PICC™	
384141	22 G	2.8F	50 cm	S/L First PICC™	Bard StatLock®

Silicone Midline Catheters (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384357	20 G	3F	20 cm	S/L First Midcath Midline	
384304	20 G	3F	20 cm	S/L First Midcath Midline	BD Introsyte™ AutoGuard Introducer
384358	18 G	4F	20 cm	S/L First Midcath Midline	
384308	18 G	4F	20 cm	S/L First Midcath Midline	BD Introsyte™ AutoGuard Introducer

Polyurethane PICCs (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384516	28 G	1.2F	25 cm	S/L L-Cath™ PICC	U-Wing Splittable Introducer Needle
384518	28 G	1.2F	25 cm	S/L L-Cath™ PICC Complete Kit	U-Wing Splittable Needle Introducer Needle and 384835
384539	26 G	1.9F	30 cm	S/L L-Cath™ PICC	
384559	26 G	1.9F	30 cm	S/L L-Cath™ PICC Complete Kit	BD Introsyte -N Introducer and 384835
384548	26 G	1.9F	30 cm	S/L L-Cath™ PICC Complete Kit	BD Introsyte AutoGuard Introducer and 384835
384466	20 G	2.6F	60 cm	D/L L-Cath™ PICC	
384463	20 G	3F	60 cm	S/L L-Cath™ PICC	
384464	18 G	3.5F	60 cm	S/L L-Cath™ PICC	
384467	18 G	3.5F	60 cm	D/L L-Cath™ PICC	

Polyurethane Midline Catheters (Box of 5)					
CATALOG #	GAUGE	FRENCH	LENGTH	DESCRIPTION	ADDITIONAL COMPONENTS
384620	20 G	3F	20 cm	S/L L-Cath Midline	
384631	20 G	3F	20 cm	S/L L-Cath Midline	BD Introsyte AutoGuard Introducer
384621	18 G	3.5F	20 cm	S/L L-Cath Midline	
384635	18 G	3.5Fr	20 cm	S/L L-Cath Midline	BD Introsyte AutoGuard Introducer

Stabilization Devices (Box of 50)	
CATALOG #	DESCRIPTION
384807	Hub Guard Foam Strips

Argon Clinical Education

Setting the standards of excellence

In today's challenging clinical environment, educational support is paramount. Argon offers neonatal, pediatric, and adult PICC placement training programs.

Argon Educational Offerings Include:

- CE credits offered for all classes
- Classes offered at your facility
- Extended Dwell Catheters for Adult and Pediatric Patients
- Extended Dwell Catheters for the Neonate
- Advanced Workshops to Expand your Competencies
- Care and maintenance of Extended Dwell Catheters

Argon Medical Devices Extended Dwell Catheter clinical education program is accredited by the California Board of Registered Nursing.

For more information please contact your Argon EDC Sales Representative.



1846 5th Concession W.,
P.O. Box 7, Rockton, ON LOR 1X0
Tel.: 519-622-4030
Toll: 1-800-996-6674
Fax: 519-622-142
mmsi@mcarthurmedical.com
www.mcarthurmedical.com