

EyeMax 2

SENSING

ANALYSIS

DELIVERY

*"Optimum comfort and protection
during phototherapy treatments."*




EYEMAX 2 OVERVIEW

Our latex-free **EyeMax 2** provides optimum comfort and protection for patients during phototherapy treatments. The one-piece, wraparound design has two attachment points that can be adjusted independently; this prevents unwanted movement and allows for a perfect fit. The ergonomic eye-protection pad was specially designed to prevent ultraviolet light leakage.

PART NUMBERS

EyeMax 2 Regular (20 pack) R300P01
EyeMax 2 Premie (20 pack) R300P02
EyeMax 2 Micro (20 pack) R300P03



1. Velcro Strapping

Adjustable band for optimum comfort

2. Eye Pads

Bigger, less pressure, and better UV protection

3. Latex Free

Helps to prevent irritant contact dermatitis



REGULAR (12.6"-14.9")
(Occipital Frontal Circumference)



PREEMIE (10.4"-12.6")
(Occipital Frontal Circumference)



MICRO (7.87"-10.4")
(Occipital Frontal Circumference)

Specifications

ML# 182 REV. D

Regular	12.6"-14.9" (Occipital Frontal Circumference)
Preemie	10.4"-12.6" (Occipital Frontal Circumference)
Micro	7.87"-10.4" (Occipital Frontal Circumference)



1846 5th Concession W., P.O. Box 7, Rockton, Ontario LOR 1X0
 Tel.: (519) 622-4030 • Toll Free: 1-800-996-6674 • Fax: (519) 622 - 1142
 Email: mmsi@mcarthurmedical.com • Web: mcarthurmedical.com