




Regular I.V. Pole



Pediatric I.V. Pole



Folding I.V. Pole  
With Bag

## I.V. POLE SOLUTIONS

The Provita® Regular I.V. Poles from Maxtec allow a low-cost alternative for patient care stands.

The Pediatric pole combines a colorful presence with a creative design to help create a more comfortable setting in the patient-area.

The Folding stand is ideal for anywhere space-savings during storage is a concern, or for Homecare applications.

## PART NUMBERS

- Normal I.V. Pole . . . . . R100P85-001
- Pediatric I.V. Pole . . . . . R100P85-002
- Folding I.V. Pole . . . . . R100P85-003

---

## Normal, Specifications


Material:	Chromed steel
Color of plastic parts:	Blue
Bottle holder:	Plastic 4 safety bottle hooks - 2 kg / hook
Adjustment:	Screw safety height adjustment
Regulating range:	1.350 - 2.150 mm
Tube diameter:	25 / 18 mm
Length of outer tube:	1.000 mm
Base diameter:	635 mm
Base leg profile:	25 x 25 mm
Casters:	Twin casters Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable
Total- / base weight (net):	Appr. 3,5 kg / 1,9 kg



---

## Pediatric, Specifications

Material:	Aluminum / plastic
Color of plastic parts:	Yellow / blue / green
Bottle holder:	Plastic 4 safety bottle hooks - 2 kg / hook
Adjustment:	Screw safety height adjustment
Regulating range:	1.350 - 2.000 mm
Tube diameter:	25 / 18 mm
Base diameter:	635 mm
Casters:	Twin casters Ø 50 mm with airplane casters
Total / base weight (net):	Appr. 3,0 / 1,3 kg



---

## Folding, Specifications

Material:	Aluminum / plastic
Bottle holder:	2 safety bottle hooks - 0,5 kg / hook
Adjustment:	Screw safety height adjustment
Regulating range:	830 - 1.780 mm
Tube diameter:	20 / 15 / 11 mm
Base diameter:	680 mm
Total load capacity:	3 kg
Total weight (net):	appr. 2 kg
Features:	Anti-slide stopper base, including carrying bag



ML# 280 REV. A