

home  
health  
scales



model RL-440HH



model RL-330HHD

model RL-330HHL

**RICE LAKE**  
HEALTH SOLUTIONS

To be the best by every measure<sup>®</sup>


[www.ricelake.com/health](http://www.ricelake.com/health)

## home health scales

RICE LAKE'S HOME HEALTH SCALES are ideal for monitoring and maintaining your weight at home. Choose from our professional digital or dial models to meet your needs.

The **RL-440HH** features the latest digital technology in a sleek professional form. The 1.5 inch, bright LCD displays weight in pounds, kilograms or stones. The auto hold function continues to show weight value on the display, even after stepping off of the scale. Tap to turn on; it automatically switches off. Lithium battery is included.

The **RL-330 Series** includes the HHL that displays weight in pounds, and the HHD that displays in pounds and kilograms. Each steel base features a large non-slip platform and a nearly 8-inch diameter dial. The rounded form can be positioned in a corner providing easy access. Bright markers can be used for record and goal keeping.



model RL-330HHL

model RL-330HHD

### RL-440HH

#### Standard Features

- Dimensions: 13.7 in x 12.9 in x 2 in
- 440 lb x 0.2 lb (200 kg x 0.1 kg)  
1 st - 10 st x 0.5 lb; 10 st - 31 st x 1.0 lb (*st = stone*)
- Display weight in pounds, kilograms or stones
- Easy to read 1.5 inch LCD display
- Durable platform with non-slip surface
- Auto hold function
- Lithium battery included
- Tap to turn on, automatic switch off

#### Warranty

- One year limited warranty

### RL-330 Series

#### Standard Features

- Dimensions: 16.93 in x 11.22 in x 3.3 in
- Easy to read 7 7/8 inch diameter dial
- Colorful markers for record purpose
- Durable steel base and platform
- Non-slip mat

#### Model RL-330HHL

- Weighs in pounds
- 330 lb x 1 lb

#### Model RL-330HHD

- Weighs in pounds or kilograms
- 330 lb (150 kg) x 1 lb(0.5 kg)

#### Warranty

- One year limited warranty

	RL-440HH	RL-330HHL	RL-330HHD
<b>Capacity</b>	440 lb (200 kg) (31 st)	330 lb	330 lb (150 kg)
<b>Graduations</b>	0.2 lb (0.1 kg) 1 st - 10 st x 0.5 lb; 10 st - 31 st x 1.0 lb	1 lb	1 lb (0.5 kg)

Your Rice Lake Weighing Systems distributor is:

**Canadian Distributor:**

McArthur Medical Sales Inc.  
1846 5th Concession West  
Rockton, ON L0R 1X0  
T+ 800.996.6674 519.622.4030 F+ 519.622.1142  
mmsi@mcarthurmedical.com  
www.mcarthurmedical.com

