

King Vision®

Video laryngoscope



Distributed by
Ambu

Key Benefits

Portable

- Disposable blades allow economical use of the King Vision video laryngoscope for all intubations.
- With the King Vision you can improve patient care by having high performance visualization capabilities for every intubation, at an affordable price.

Affordable

- The King Vision is light weight, self-contained, battery operated and water resistant.
- The King Vision can be used continuously for more than 1 hour making it highly portable.
- Blades are individually packaged. Each blade is disposable, eliminating concerns on cross contamination.

Durable

- Designed to be your primary tool for all intubations.
- Video output capability is compatible with external monitor and recording devices.
- The robust, full-color, non-glare display can resist repeated cleaning and normal use wear and tear.

King Vision® video laryngoscope

Visionary performance in the palm of your hand. We know that reliable performance is a necessity, not an option. When it comes to visualising the airway for routine and difficult intubations, you need a tool that you can depend on each and every time. The King Vision video laryngoscope is designed with your expectations in mind.

Portable


The King Vision display is light weight, self-contained, battery operated, and water resistant. The reusable display comes packaged in a protective, foam case. Blades are individually packaged making the King Vision very portable.

Affordable

The disposable blades allow economical use of the King Vision for all of your intubations. The low cost per procedure combined with the high performance visualization capabilities make this the perfect value for the best patient care. Take the guess-work out of your intubations. See the vocal cords for every placement. Use the King Vision for all of your intubations.

Durable

The King Vision is designed to be your primary tool for intubations. The high-quality, made in America, display comes with a 1-year warranty. The robust, full-color, smudge-proof, non-glare display can resist repeated cleaning and normal use wear and tear.



REUSABLE DISPLAY



DISPOSABLE CHanneled BLADE



DISPOSABLE STANDARD BLADE

Ordering specifications

Item no.	Type	Quantity
KVLO3	Standard Blade, Size 3	10
KVLO3C	Channeled Blade, Size 3	10
KVISO1	Reusable Display	1
KVLKIT3	Vision Display, Case, AAA Batteries, 3 Channeled and 1 Standard Size 3 Vision Blades	1
KVCABL	Custom Video-out Cable (2.75 m / 9')	1

Ambu USA
6230 Old Dobbin Lane
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
www.ambuUSA.com

CE 0086 **US: Rx only**