

OPTIONAL ACCESSORIES



## 30 WATT SOLAR PANEL

Portable, Foldable, Rugged Charging

### PART NUMBERS

Mfg. #	Description
EXP48PRO	LiFePO4 Battery System 48,000mAh
EXP96PRO	LiFePO4 Battery System 96,000mAh
EXP-SC30W	30W Portable Solar Panel
S1-60	System One REMstar 60 Series DC Power Cord
DS-12V	Respironics DreamStation DC Power Cord
RM-9	Resmed S9 DC Power Cord
RM-10	Resmed S10 DC Power Cord

Compatible with most CPAP devices including  
**Philips Respironics and ResMed**

# EXP PRO

Manufactured by Western Medical Inc.  
Offered by:

**MMSI**  
McARTHUR  
MEDICAL SALES INC.

1846 5th Concession W., P.O. Box 7, Rockton, Ontario LOR 1X0  
Tel.: (519) 622-4030 • Toll Free: 1-800-996-6674 • Fax: (519) 622 - 1142  
Email: [mmsi@mcarthurmedical.com](mailto:mmsi@mcarthurmedical.com) • Web: [mcarthurmedical.com](http://mcarthurmedical.com)

# EXP PRO

POWER TO EXPLORE



## PORTABLE HANDHELD POWER

POWER WHEREVER YOU ARE

## LIFEPO4

# 10X

MORE LIFE  
CYCLES

# SLEEP THERAPY WHEREVER YOU ARE

The EXP PRO Batteries power your CPAP machine anywhere at anytime. They're lightweight, portable and can power your machine and other devices for days on end. You'll never lose sleep again.

Camping  
Power Outages  
Travel Abroad



## EXP48PRO

**2-3  
NIGHTS**  
CPAP Usage



- TSA COMPLIANT**
- LED Light
- Type A&C +QC
- Runs **1** CPAP Machine\*

# PREMIUM LITHIUM TECHNOLOGY: LiFePo4

Not all lithium batteries are the same. EXP PRO Battery systems are engineered with LiFePo4 chemistry, an advanced lithium technology that optimizes power and performance so you can sleep easy, wherever you are.\*

10X Life Cycles ≤ 5,000+  
2 to 3X Faster Charging  
6X Discharge Power  
65% Lighter Weight  
Extreme Weather Performance

## EXP96PRO

**4-6  
NIGHTS**  
CPAP Usage



- LED Light
- Type A&C +QC
- Runs **2** CPAP Machines\*

*\* performance varies by model power draw*

	EXP48PRO	EXP96PRO
CPAP Nights	2-3	4-6
CPAP Machines	1	1-2
Capacity	154Wh/12.8Ah @ 12V	307Wh/25.6Ah @ 12V
Power Ports	15A Max	15A Max
USB Ports	2 Type-A QC/1 Type-C	2 Type-A QC/1 Type-C
Solar Input	18-25V @ 2.0 Amps	18-25V @ 2.0 Amps
Charge Time	4-6 Hrs	8-12 Hrs
Weight	3.25lbs	6.5 lbs
HWL	1.6" x 6" x 8.2"	3" x 6" x 8.2"
Warranty	2 Years	2 Years



\* compared to lead acid batteries