

FIBER OPTIC GREENLINE® / D™ BLADE

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

SunMed's Fiber Optic GreenLine®/D™ disposable laryngoscope blades help solve contamination problems and eliminate the cost and time spent cleaning blades and returning them for use. These stainless steel blades also answer the professional's request for a non-plastic disposable blade. The polished acrylic stem produces exceptional illumination. The extended safety heel prevents blade-to-handle contact. All sizes and styles are compatible with Green System fiber optic handles and are ISO 7376 compliant.



MILLER GREENLINE® / D™

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5332-01	Infant	1	92 mm	22 mm	20
5-5332-02	Child	2	100 mm	21 mm	20
5-5332-03	Medium Adult	3	130 mm	22 mm	20
5-5332-35	Ext. Med. Adult	3.5	145 mm	22 mm	20
5-5332-04	Large Adult	4	155 mm	25 mm	20

MACINTOSH GREENLINE® / D™

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5333-20	Premature	00	67 mm	10 mm	20
5-5333-00	Newborn	0	77 mm	10 mm	20
5-5333-01	Infant	1	102 mm	10 mm	20
5-5333-02	Child	2	153 mm	12 mm	20
5-5333-03	Medium Adult	3	195 mm	14 mm	20
5-5333-04	Large Adult	4	205 mm	18 mm	20

