

What is OxyGo?

The OxyGo® family of portable oxygen concentrators (POC) is designed to keep oxygen patients—like you—GOING!

An OxyGo POC is a lightweight medical device that takes the air around you and converts it into the oxygen prescribed by your doctor.

The OxyGo family has three members:

OxyGo NEXT™
OxyGo and OxyGo FIT™

All are industry leaders in reliability, usability and quality. And all are clinically validated to be used 24/7... so they can supply the medical oxygen you need at home or on the go.

OxyGo is different from other POCs because it shows you exactly how much time you have left on the battery, and you can change the battery easily. You can also charge OxyGo and use at the same time—either in the car or at home!

Will OxyGo work for me?

OxyGo NEXT, OxyGo and OxyGo FIT are available by prescription only. If you have been prescribed oxygen, there is a good chance that a member of the OxyGo family will work for you. However, you will need your oxygen provider to perform an oxygen assessment to see if you can saturate with OxyGo.

**OxyGo® is ready to go
when you are!**

Whether you are out on the town or home with family—now you can always have your oxygen with you. Learn more at: www.oxygo.life

Your OxyGo Provider:

Are you a medical professional?

OxyGo has been clinically tested to saturate patients 24/7... and its pulse dose is sensitive enough to provide most patients with the clinically effective prescription required.

OxyGo is available with local support through your area oxygen provider.

For more information: www.oxygo.life/clinical

Follow Us!



**MADE
IN THE USA**

Battery life statistics based on average values may vary up to +/- 10% depending on patient use.

PN B117-0819

www.oxygo.life

Now with the **NEW** “My OxyGo” APP!



You Need Oxygen... NOT TANKS.



OxyGO Keep Going.™

www.oxygo.life



NEW!
My OxyGo
App*



OxyGo NEXT™

7.19" x 3.26" x 8.15" 4.7 lbs.

OxyGo NEXT is the next generation of portable oxygen. Its ground-breaking size and battery life have never been seen before at this high an oxygen dose, and it is setting the standard for portable oxygen concentrators. OxyGo NEXT has an easy to read LCD display, simple controls, and can be charged at home or on the go—giving you the freedom to embrace what's NEXT!

OxyGo NEXT

(with 8-cell battery) #1400-3000

- 4.7 lbs. | 7.19" x 3.26" x 8.15"
- 6 Pulse Flow Settings (1260 ml/min)
- Up to 13 hour battery life*
- Only 38 dB. at setting 2
- Included 5-year warranty on the unit, 1-year warranty on sieve beds & battery
- Bluetooth-enabled unit connects to your mobile device with the **My OxyGo App**.** Check your battery life, settings and more from your phone!

OxyGo NEXT Comes With:

- Standard, 8-cell battery
- DC car charger, and AC wall charger
- Nasal cannula
- Carry bag with padded, adjustable strap



NEW!
My OxyGo
App*



OxyGo FIT™

5.9" x 2.7" x 7.2" 2.8 lbs.

OxyGo FIT weighs just 2.8 pounds and has 3 settings, with 210 ml of oxygen per setting. The OxyGo FIT is well constructed and will operate for almost 3 hours using a standard, user-replaceable battery. It also has an optional larger battery, which will power the unit for up to 5 hours* and adds less than one-half pound in weight.

OxyGo FIT

(with 4-cell battery) #1400-2000

- 2.8 lbs. | 5.9" x 2.7" x 7.2"
- 3 Pulse Flow Settings (630 ml/min)
- Up to 5 hour battery life*
- Only 40 dB. at setting 2
- Included 3-year warranty on the unit, 1-year warranty on sieve beds & battery
- Bluetooth-enabled unit connects to your mobile device with the **My OxyGo App**.** Check your battery life, settings and more from your phone!

OxyGo FIT Comes With:

- Standard, 4-cell battery
- DC car charger, and AC wall charger
- Nasal cannula
- Carry bag with padded, adjustable strap



OxyGo 5 Setting

8.75" x 3.0" x 7.25" 4.8 lbs.

The original OxyGo has 5 pulse flow settings, 1,050 ml output, and up to 10.25 hours* of battery life—all in a compact, portable unit that weighs less than 5 pounds. OxyGo also has user-replaceable batteries and sieve beds, so it's ready when you need it. And with its 5 year warranty, OxyGo is truly designed to help you Keep Going.

OxyGo 5 Setting

(with 8-cell battery) #1400-1000

- 4.8 lbs. | 8.75" x 3.0" x 7.25"
- 5 Pulse Flow Settings (1050 ml/min)
- Up to 10.25 hour battery life*
- Only 39 dB. at setting 2
- Included 5-year warranty on the unit, 1-year warranty on sieve beds & battery

OxyGo 5 Setting Comes With:

- Standard, 8-cell battery
- DC car charger, and AC wall charger
- Nasal cannula
- Carry bag with padded, adjustable strap

* double cell battery, setting 1

** Available on OxyGo FIT POCs produced after 2/19 and OxyGo NEXT