

Pulmodyne®



GO-PAP™

with Integrated Nebulization

- Run time of approximately **40 mins** on a D-Cylinder Tank
- **PEEP settings of 5/7.5/10**
- Receive up to **90 LPM** to your CHF or COPD Patients who need it most
- **Set it at 10LPM and let it flow...**

MMSI
McARTHUR
MEDICAL SALES INC.

1846 5th Concession W.,
P.O. Box 71 Rockton, ON LOR 1X0
Tel.: 519-622-4030
Toll: 1-800-996-6674
Fax: 519-622-142
mmsi@mcarthurmedical.com
www.mcarthurmedical.com