

# Immersion Swaddle Bathing

Transform bathing to a positive experience



## Benefits of Immersion Swaddle bathing:

- Minimize temperature loss <sup>1-3</sup>
- Decreases crying <sup>1,3</sup>
- Reduces physical and behavioral stress <sup>2, 4-6</sup>
- Supports family centered care <sup>6,7,10</sup>
- Enhances ability to feed after bath <sup>6,10</sup>
- Evidence-Based Practice <sup>1-6, 8-10</sup>

Swaddle bathing with the TurtleTub complements delayed bathing and minimizes adverse clinical outcomes associated with traditional bathing methods. After participating in hospital bathing, parents can continue to swaddle bathe at home. Preterm infants, full-term infants, and infants with NAS all benefit from immersion swaddle bathing. In addition, the TurtleTub is a strong parent pleaser, providing a marketing advantage for your hospital.



**CATAPULT**  
PRODUCTS

[www.catapult-products.com](http://www.catapult-products.com)

Early Experiences Matter

# The TurtleTub System:

Designed for infant safety, comfort, ease of use and infection control.



## TurtleTub Infant Bathtubs:

- Ideal for swaddle bathing and immersion tub bathing
- Options for single-patient and multi-patient use
- Ribs and high sides provide excellent support
- Integrated temperature strip
- Smooth contours for easy cleaning
- Convenient pour spot
- Stackable
- Complies with ASTM F2670 standard for infant bathtubs

## TurtleTub Dimensions:

22"L x 18" W x 6"H

## TurtleTub Disposable Bath Liners (for multi-patient tub use):

- Strong/thick poly with elastic holds liner securely in TurtleTub
- Infection control compliant
- Non-leaching poly material complies to USP Physiochemical 661 and Cytotoxicity tests providing a safe bath
- Available in packs with accessories or packaged individually



## TurtleTub Swaddle Bathing Blankets:

- Reusable fleece blankets with excellent thermal properties
- Disposable non-woven for hospitals without laundry options
- Reusable cotton blankets provide a cost-effective fleece alternative

## TurtleTub Bathing Cart:

- 4-sided guard rail keeps TurtleTub secure
- Two locking wheels for safe bathing
- Two smooth gliding “soft-close” drawers

### Specifications:

- 20” wide x 16” deep x 34” height (not including guard rail)
- 4-sided guard rail (adds 3.5” to height)
- Stainless steel grade 304

Cart Caddy hangs conveniently on the guard rail to hold bathing supplies



Optional Cart Caddy  
(sold separately)

## TurtleTub Accessories:

Available individually and in our kits and packs  
(See ordering info for details)

### Baby Wash:

- Single use packets (10mL)
- pH neutral and mild
- Made without dyes, sulfates, parabens or phthalates
- Fragrance-free



22 oz. rinse cup



Mesh storage bag to hang, dry  
and store the TurtleTub



Brush and comb set for  
full-term infants



Small scalp brush for  
small or premature infants

## TURTLETUB™ ORDERING INFORMATION

Item #	Description	Qty/Case
<b>Complete Swaddle Bathing Kits</b>		
<b>TT1100</b>	<b>TurtleTub Standard Kit</b> Each kit contains: 1 TurtleTub with Integrated Temperature Indicator; 1 Neutral PH Mild Baby Wash; 1 Small Fleece Blanket; 1 Rinse Cup; 1 Scalp brush; 6 Washcloths & Instructions for Use.	10 Kits
<b>TT1100MB</b>	<b>TurtleTub Mother Baby Unit Standard Kit</b> Each kit contains: 1 TurtleTub with Integrated Temperature Indicator; 1 Neutral PH Mild Baby Wash; 1 Large Cotton Flannel Blanket; 1 Rinse Cup; 1 Brush & Comb set; 6 Washcloths & Instructions for Use.	10 Kits
<b>TT1106</b>	<b>TurtleTub NICU Starter Kit w/Liner</b> - Each kit contains: 5 TurtleTub Bathtubs (each with Integrated Temperature Indicator & Rinse Cup); 10 Small Fleece Blankets and 25 Accessory Packs w/Liner (TT1107).	1 Kit
<b>TT1113</b>	<b>TurtleTub Deluxe Kit</b> - Each kit contains: 1 TurtleTub with Integrated Temperature Indicator; 1 Neutral PH Mild Baby Wash; 1 Large Fleece Blanket; 1 Rinse Cup; 1 Scalp brush; 1 Mesh Storage Bag; 6 Washcloths & Instructions for Use.	10 Kits
<b>TurtleTub - Bathtubs (Tub only - no accessories)</b>		
<b>TT1104</b>	<b>TurtleTub Bathtub</b> (with Integrated Temperature Indicator)	10
<b>TT1104-2</b>	<b>TurtleTub Bathtub</b> (with Integrated Temperature Indicator)	2
<b>Disposable Bath Liners &amp; Liner Packs</b>		
<b>TT1124</b>	<b>Mother Baby Unit Premium Bath Accessory Pack</b> Each pack contains: 1 Disposable Liner; 1 Disposable Swaddle Bath Blanket; 1 Neutral pH Mild Baby Wash; 1 Brush & Comb set & 6 Washcloths.	50
<b>TT1125</b>	<b>Mother Baby Unit Deluxe Bath Accessory Pack</b> Each pack contains: 1 Disposable Liner; 1 Disposable Swaddle Bath Blanket; 1 Brush & Comb set & 6 Washcloths.	50
<b>TT1122</b>	<b>Mother Baby Unit Standard Bath Accessory</b> Each pack contains: 1 Disposable Liner; 1 Neutral pH Mild Baby Wash; 1 Brush & Comb set & 6 Washcloths.	50
<b>TT1123</b>	<b>Mother Baby Unit Basic Bath Accessory Pack</b> Each pack contains: 1 Disposable Liner; 1 Brush & Comb set & 6 Washcloths.	50
<b>TT1107</b>	<b>NICU Standard Bath Accessory Pack</b> - Each pack contains: 1 Disposable Liner; 1 Neutral pH Mild Baby Wash; 1 Scalp brush & 6 Washcloths.	50
<b>TT1126</b>	<b>NICU Basic Bath Accessory Pack</b> Each pack contains: 1 Disposable Liner; 1 Scalp brush & 6 Washcloths.	50
<b>TT1105</b>	<b>TurtleTub Disposable Bath Liner</b> - Individually Packed	100
<b>Swaddle Bathing Blankets</b>		
<b>TT1103</b>	<b>Fleece Reusable Swaddle Bathing Blanket - Small</b> (20" x 23")	25
<b>TT1110</b>	<b>Fleece Reusable Swaddle Bathing Blanket - Large</b> (20" x 30")	25
<b>TT1115</b>	<b>Cotton Flannel Reusable Swaddle Bathing Blanket - Large</b> (20" x 30")	25
<b>TT1119</b>	<b>Non-woven Disposable Swaddle Bathing Blanket - Large</b> (20" x 30")	100
<b>Mobile Bathing Cart</b>		
<b>TT1114</b>	<b>Stainless Steel 2-drawer Bathing Cart</b> - 20" W x 16" D x 34" H, 4-sided guard rail adds 3.5" to Height	1
<b>TT1120</b>	<b>Cart Caddy</b> - 14.5" x 3.75" x 3.75"	1
<b>Baby Wash, Brush &amp; Comb Set, Storage Bag &amp; Rinse Cup</b>		
<b>TT1102</b>	<b>Baby Wash Single Use Packets</b> - 10mL Neutral pH Mild Baby Wash	100
<b>TT1127</b>	<b>Brush &amp; Comb Set</b>	100
<b>TT1108</b>	<b>TurtleTub Storage Bags</b>	10
<b>TT1109</b>	<b>TurtleTub Rinse Cups</b> - 22 oz	25

### REFERENCES

- Edraki M, et al. Comparing the effects of swaddled and conventional bathing methods on body temperature and crying duration in premature infants: a randomized clinical trial. *Journal of Caring Sciences*. 2014; 3 (2): 83-91.
- Caka SY, Gözen D. Effects of swaddled and traditional tub bathing methods on crying and physiological responses of newborns. *J Spec Pediatr Nurs*. 2018;23:e12202
- Swapna, G., et al. Relative Effectiveness of Swaddle Bath and Conventional Bath on Level of Thermal Stability and Crying Duration among Preterm Infants at Selected Hospital in North India. *ICCRJNR*. Jan-Jun 2017; 2(1): 34-54.
- Ceylan SS, Bolisik. Effects of swaddled and sponge bathing methods on signs of stress and pain in premature newborns: Implications for Evidence-Based Practice. *World Views on Evidence Based Nursing*. 2018; 15: 296-303. doi 10.1111/wvn.12299
- Paran M, Edraki M, Montaseti S, Razavi Nejad M. Comparing the Effects of Swaddle and Conventional Bathing Methods on Behavioral Responses in Preterm Neonates *IJN*. 2016; NOV, 7(4) DOI:10.22038/ijn.2016.7778
- Hall K. Practising developmentally supportive care during infant bathing: reducing stress through swaddle bathing. *Infant*. 2008; 4(6): 198-201.
- Craig, JW, et al. Recommendations for involving the family in developmental care of the NICU baby. *Journal of Perinatology*. 2015; 35(Suppl 1): S5-S8.
- Quirashy K, et al. A protocol for swaddled bathing in the neonatal intensive care unit. *NAINR*. 2013; (1): 48-50.
- Association of Women's Health, Obstetric, and Neonatal Nurses (AWHONN). *Neonatal Skin Care: Evidence-based Clinical Practice Guideline*. 2018; 41-53.
- Fern D, et al. Swaddled bathing in the newborn intensive care unit. *Newborn and Infant Nursing Reviews*. 2002; 2(1), 3-4.

©2021 Catapult Products LLC. TurtleTub is a Trademark of Catapult Products LLC. TurtleTub is protected by US and foreign patents pending.



Manufactured by:  
**CATAPULT PRODUCTS, LLC**  
 4677 S. Cherry St.  
 Salt Lake City, UT 84123

Distributed by:



**1846 5th Concession W., P.O. Box 7, Rockton, Ontario LOR 1X0**  
**Tel.: (519) 622-4030 • Toll Free: 1-800-996-6674 • Fax: (519) 622 - 1142**  
**Email: mmsi@mcarthurmedical.com • Web: mcarthurmedical.com**