

# Zero-Resistance™ Options+™ Bottle System Specialty Feeding System Starter Kit

Designed to assist in the management and treatment of complicated oral feeding issues requiring a compression nipple. This may include conditions including cleft lip/palate, tongue tie, lip tie, oro-neuromotor dysfunctions, congenital heart disease, or craniofacial anomalies.

- Promotes “Infant-Directed” feedings
- Typical feeding bottle design; fully reusable
- Helps decrease gassiness/air intake and burping with non-squeeze method
- Assembled with an Infant-Paced Feeding Valve and Level 1 Narrow Nipple
- Compatible with any of the Dr. Brown’s® eight (8) Narrow Nipples
- Individually wrapped
- BPA free
- Infant-Paced Feeding Valves are not available for individual purchase
- Each bottle system purchased comes with 2 Infant-Paced Feeding Valves
- The full bottle system is dishwasher, microwave steam bag, and steam sterilization (autoclave) safe
- Preferred by parents and medical professionals

\*Autoclave recommended not to exceed 30 cycles and is NOT recommended for the Infant-Paced Feeding Valve.

## RETAIL PACKAGE

ITEM #	ITEM DESCRIPTION	RETAIL UPC	# PER CASE
SB05007-WEB	Zero-Resistance™ Specialty Feeding System Starter Kit	0-72239-30727-4	6

MASTER CARTON: 6 UNITS/CASE, Frustration-Free packaging

**Set Includes:** 3 - 8oz/250 mL bottles, 2 - 4 oz/120 mL bottles, 1 bag of 5 additional valves, 3 - 5.5" wire brushes, 2 storage caps, 1 bag of 5 travel disks, 1 bag of 5 Level 2 Nipples

	DIMENSIONS			WEIGHT			
	WIDTH	LENGTH	HEIGHT	Cubic Ft	10 <sup>th</sup> 's	POUNDS	OUNCES
BULK	14.9"	15.4"	10.3"	1.37 cu ft	8.32 lb	8.0 lb	5.1 oz

Included Documents: Specialty Feeding System Consumer Instruction Insert



**For optimal efficacy of this product, the Specialty Feeding System must be assembled with all 3 components –**

- **Dr. Brown’s® Nipples – Quality controlled to assure dependable, reliable flow rates**
- **Infant-Paced Feeding Valve**
- **Internal Vent System**

1846 5th Concession W., P.O. Box 7, Rockton, Ontario LOR 1X0  
Tel: 519.622.4030 • Toll Free: 1.800.996.6674 674 • Fax: 519.622.1142  
Email: mmsi@mcARTHURmedical.com • Web: mcARTHURmedical.com

