


# ODAPT<sup>®</sup>

O<sub>2</sub>

An innovative solution designed to seamlessly integrate an oxygen analyzer with the most commonly used sizes of oxygen tubing.



Swivel connection

Analyzer port

22mm OD

15 mm OD

1/4 Cannula Tubing

# ODAPT<sup>®</sup>

O2

Lightweight

Versatile

Multifunction



## ORDERING INFORMATION

ITEM#	DESCRIPTION	QTY
420-1980	ODAPT O2	25/box

ODAPT is a registered trademark used under license by McArthur Medical Sales Inc. US and International Patents Granted and Pending.

McArthur Medical Sales Inc  
 1846 5th Concession Rd W  
 Rockton, ON L0R 1X0  
 1-800-996-6674  
 marketing@mcarthurmedical.com  
 www.mcarthurmedical.com



**MMSI**  
**McARTHUR**  
 MEDICAL SALES INC.

REV.1 ODA20240206