

A newborn baby is lying in a hospital bed, wearing a nasal cannula. The baby's face is the central focus, with the cannula clearly visible. The background is a blurred hospital room. The text 'NUBRE' is overlaid at the top in a large, blue, outlined font, with a trademark symbol. Below it, 'by MCARTHUR MEDICAL' is written in a smaller, blue, sans-serif font. The main title 'Enteral Feeding Solutions' is centered in a large, bold, blue font. At the bottom, a descriptive sentence is written in a smaller, bold, blue font.

NUBRE™

by MCARTHUR MEDICAL

Enteral Feeding Solutions

NUBRE™ is an enteral feeding solution that offers a blend of nutrients tailored to support patients nutritional needs with precision and efficiency

Feeding Tubes

- NUBRE™ feeding tubes are available with ENFit® connectors
- Feeding tube sizes range from 4FR to 12FR
- Available in Orange and Purple to suit all customer's needs
- Available in Polyurethane (PU) and PVC materials
- Free of DEHP, BPA and Latex
- Feeding tubes have a radiopaque strip and bold markings
- Sterile

ENFit® Sterile PU Feeding Tube

Orange	Purple	Description	UOM
MM1404040	MM1404040P	4Fr 40cm	BX/10
MM1405040	MM1405040P	5Fr 40cm	BX/10
MM1406040	MM1406040P	6Fr 40cm	BX/10
MM1408040	MM1408040P	8Fr 40cm	BX/10
MM1404050	MM1404050P	4Fr 50cm	BX/10
MM1405050	MM1405050P	5Fr 50cm	BX/10
MM1406050	MM1406050P	6Fr 50cm	BX/10
MM1408050	MM1408050P	8Fr 50cm	BX/10
MM1404080	MM1404080P	4Fr 80cm	BX/10
MM1405080	MM1405080P	5Fr 80cm	BX/10
MM1406080	MM1406080P	6Fr 80cm	BX/10
MM1408080	MM1408080P	8Fr 80cm	BX/10
MM14010080	MM14010080P	10Fr 80cm	BX/10
MM14012080	MM14012080P	12Fr 80cm	BX/10
MM1404090	MM1404090P	4Fr 90cm	BX/10
MM1405090	MM1405090P	5Fr 90cm	BX/10
MM1406090	MM1406090P	6Fr 90cm	BX/10
MM1408090	MM1408090P	8Fr 90cm	BX/10
MM140100110	MM140100110P	10Fr 110cm	BX/10
MM140120110	MM140120110P	12Fr 110cm	BX/10

ENFit® Sterile PVC Feeding Tube

Orange	Purple	Description	UOM
MM1444040	MM1444040P	4Fr 40cm	BX/25
MM1445040	MM1445040P	5Fr 40cm	BX/25
MM1446040	MM1446040P	6Fr 40cm	BX/25
MM1448040	MM1448040P	8Fr 40cm	BX/25
MM1445050	MM1445050P	5Fr 50cm	BX/25
MM1446050	MM1446050P	6Fr 50cm	BX/25
MM1448050	MM1448050P	8Fr 50cm	BX/25
MM1444080	MM1444080P	4Fr 80cm	BX/25
MM1445080	MM1445080P	5Fr 80cm	BX/25
MM1446080	MM1446080P	6Fr 80cm	BX/25
MM1448080	MM1448080P	8Fr 80cm	BX/25
MM14410080	MM14410080P	10Fr 80cm	BX/10
MM14412080	MM14412080P	12Fr 80cm	BX/10
MM1445090	MM1445090P	5Fr 90cm	BX/25
MM1446090	MM1446090P	6Fr 90cm	BX/25
MM1448090	MM1448090P	8Fr 90cm	BX/25
MM144100110	MM144100110P	10Fr 110cm	BX/25
MM144120110	MM144120110P	12Fr 110cm	BX/25

The NUBRE™ feeding tube is intended for nasogastric or orogastric delivery of various types of liquid nutritional media through the gastrointestinal tract of neonatal and small pediatric patients. The Polyurethane feeding tube is intended for use up to 30 days. The PVC feeding tube is intended for use up to 72 hours.



Tel. (519) 622-4030 • Toll Free Line 1-800-996-6674 • Fax (519) 622-1142
www.mcarthurmedical.com • mmsi@mcarthurmedical.com
 Distributed by McArthur Medical Sales Inc