



Arterial Blood Sampling

Product Catalog

Canada

Table of Contents

Arterial Blood Sampling

Pro-Vent Plus™ 3 mL Arterial Blood Sampling Kits–Full.....	2
Pro-Vent Plus™ 3 mL Arterial Blood Sampling Kits–Strip.....	3
Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits–Full.....	3
Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits–Strip.....	4
Pulsator™ Arterial Blood Sampling Full and Strip Kits.....	4
Portex™ LineDraw™ Plus Arterial Blood Sampling Kits Dry Heparin.....	5

Arterial Blood Sampling

Pro-Vent Plus™ 3 mL Arterial Blood Sampling Kits– Full*

100 Units Per Case	Product Codes									
Components	4589P-1	4591P-1	4599P-1	4649P-1	4650P-1	4651P-1	4660P-1	4699P-1	4699LH	G1702
3 mL Luer Slip Syringe				•	•	•	•	•	•	
3 mL Luer Lock Syringe	•	•	•							•
Needle										
22G x 1" Pre-Attached to Syringe		•								
22G x 1" Pre-Attached with Needle-Pro™ to Syringe				•	•					
22G x 1 ½" Secondary Needle	•									
22G x 1 ½" Secondary Needle Attached to Needle-Pro			•					•	•	
23G x 1" Pre-Attached to Syringe	•									
23G x 1" Pre-Attached with Needle-Pro™ to Syringe			•			•	•	•	•	•
Filter-Pro™ Air Bubble Removal Device	•	•	•	•	•	•	•	•	•	•
Point-Lok™ Needle Protection Device	•	•								
Alcohol Prep Pad	•	•	•	•	•	•	•	•	•	•
Iodophor Prep Pad	•	•	•	•			•	•	•	
Gauze Pad (2" x 2")	•	•	•	•	•	•	•	•	•	•
Patient Label	•	•	•	•	•	•	•	•	•	•
Biohazard Labeled Ice Bag	•	•	•	•	•	•	•	•	•	•
Adhesive Bandage (1" x 3")	•	•	•	•			•	•	•	
Sterilized										
Fully Sterile	•	•	•	•	•	•	•	•	•	•

Pro-Vent Plus™ 3 mL Arterial Blood Sampling Kits–Strip*								
200 Units Per Case	Product Codes							
Components	4346-2	4588P-2	4590P-2	4598P-2	4645P-2	4646-2	4698LH	4698P-2
3 mL Luer Slip Syringe					•	•	•	•
3 mL Luer Lock Syringe	•	•	•	•				
Needle								
22G x 1" Pre-Attached to Syringe			•					
22G x 1" Pre-Attached with Needle-Pro™ to Syringe					•			
23G x 1" Pre-Attached to Syringe		•						
23G x 1" Pre-Attached with Needle-Pro™ to Syringe				•			•	•
Filter-Pro™ Air Bubble Removal Device	•	•	•	•	•	•	•	•
Point-Lok™ Needle Protection Device		•	•					
Adhesive Bandage 3/4" x 3"								•
Sterilized								
Fully Sterile	•	•	•	•	•		•	•
Sterile Fluid Path						•		

*Pro-Vent™ 3 mL syringes contain approximately 23.5 I.U. per mL of advanced formula dry lithium heparin.

*Pro-Vent™ Low Heparin (LH Codes) 3 mL syringes contain approximately 7 I.U. per mL of dry lithium heparin.

Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits–Full*	
100 Units Per Case	Product Codes
Components	4640P-1
1 mL Luer Slip Syringe	•
Needle	
25G x 5/8" Pre-Attached with Needle-Pro™ to Syringe	•
23G x 1" Secondary Needle Pre-Attached with Needle-Pro™ to Syringe	•
Filter-Pro™ Air Bubble Removal Device	•
Point-Lok™ Needle Protection Device	
Alcohol Prep Pad	•
Iodophor Prep Pad	•
Gauze Pad (2" x 2")	•
Patient Label	•
Biohazard Labeled Ice Bag	•
Adhesive Bandage (1" x 3")	•
Sterilized	
Fully Sterile	•

Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits–Strip*						
200 Units Per Case	Product Codes					
Components	4610LH	4610P-2	4616-2	4619LH	4619P-2	4629P-2
1 mL Luer Slip Syringe	•	•	•	•	•	•
Needle						
25G x 5/8" Pre-Attached to Syringe	•	•				
25G x 5/8" Pre-Attached with Needle-Pro™ to Syringe				•	•	
23G x 1" Needle Pre-Attached with Needle-Pro™ to Syringe						•
Filter-Pro™ Air Bubble Removal Device	•	•	•	•	•	•
Point-Lok™ Needle Protection Device	•	•				
Sterilized						
Fully Sterile	•	•		•	•	•
Sterile Fluid Path			•			

*Pro-Vent™ 1 mL syringes contain approximately 23.5 I.U. per mL of advanced formula dry lithium heparin.

*Pro-Vent™ Low Heparin (LH Codes) 1 mL syringes contain approximately 7 I.U. per mL of dry lithium heparin.

Pulsator™ Arterial Blood Sampling Full and Strip Kits			
	Product Codes		
Units Per Case	200	200	100
Components	4046-2	4080P-2	4084P-1
3 mL Luer Lock Syringe	•	•	•
Needle:			
22G x 1 1/2" Secondary Needle with Needle-Pro™			•
23G x 1" Pre-Attached with Needle-Pro to Syringe		•	•
Filter-Pro™ Air Bubble Removal Device	•	•	•
Alcohol Prep Pad			•
Iodophor Prep Pad			•
Gauze Pad (2" x 2")			•
Patient Label		•	•
Biohazard Labeled Ice Bag			•
Adhesive Bandage (1" x 3")		•	•
Sterilized			
Fully Sterile		•	•
Sterile Fluid Path	•		

Dry Pulsator™ 3 mL syringes contain approximately 23.5 I.U. per mL of advanced formula dry lithium heparin.

Portex™ LineDraw™ Plus Arterial Blood Sampling Kits Dry Heparin									
200 Units Per Case	Product Codes								
Components	4041-2	4041LH	4041W-2	4042-2	4042LH	4043-2	4043W-2	4040-2 Umbilical	G1777 Umbilical
1 mL Luer Slip Syringe	•	•	•						
3 mL Luer Slip Syringe						•	•		
3 mL Luer Lock Syringe				•	•			•	2x
3 mL Waste Syringe (Luer Lock)			•						
6 mL Waste Syringe (Luer Slip)							•		
Needle									
22G x 1 ½" Needle with Needle-Pro™								•	2x
Filter-Pro™ Air Bubble Removal Device	•	•	•	•	•	•	•	•	2x
Point-Lok™ Needle Protection Device								•	
Tip Cap		•	•				•		
Biohazard Labeled Ice Bag									
Sterilized									
Fully Sterile	•	•	•	•	•	•	•	•	•

Portex™ Umbilical and dry heparin Line Draw 3 mL syringes contain approximately 23.5 I.U. per mL of advanced formula dry lithium heparin. Dry heparin Line Draw Low Heparin (LH Codes) 1 mL syringes contain approximately 7 I.U. per mL of dry lithium heparin.



Pro-Vent™ Arterial Blood Sampling Kits



Portex™ LineDraw™ Plus Arterial Blood Sampling Kits



Pulsator™ Arterial Blood Sampling Kits



Filter-Pro™



Needle-Pro™



Point-Lok™

For more information or to place an order
outside of the Atlantic Provinces, please contact:

McArthur Medical Sales Inc. at 800-996-6674,
or via email at customerservice@mcarthurmedical.com
or orders@mcarthurmedical.com

For more information or to place an order
in the Atlantic Provinces, please contact:

Braemed Ltd. at 1-888-830-7744 or 902-468-7544,
or via email at info@braemed.ca.