



General Anesthesia Airway

Product Catalog

icumedical
human connections

Contents

Cuffed Magill Tubes	3
Tracheal Tube, Clear PVC, Oral, Murphy Eye, Blue Line™, Sacett™ Suction Above the Cuff, Profile™ Soft-Seal™ (Low Pressure/High Volume)	3
Tracheal Tube, Clear PVC, Murphy Eye, Oral/Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff	3
Cuffed Endotracheal tube with Murphy Eye, Clear PVC, Magill Curve.....	4
Cuffed Endotracheal tube with Preloaded Stylet, Clear PVC, Magill Curve, Murphy Eye.....	4
Tracheal Tube, Siliconized PVC, Oral/Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff.....	5
Tracheal Tube, Ivory PVC, Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff.....	5
Uncuffed Magill Tubes	6
Tracheal Tube, Siliconized PVC, Oral/Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff.....	6
Tracheal Tube, Ivory PVC, Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff.....	6
Tracheal Tube, Clear PVC, Murphy Eye, Oral/Nasal, Uncuffed	7
Tracheal Tube, Siliconized PVC, Oral/Nasal, Uncuffed	7
Tracheal Tube, Siliconized PVC, Murphy Eye, Oral/Nasal, Uncuffed.....	8
Tracheal Tube, Siliconized, Oral, Uncuffed	8
Tracheal Tube, Ivory PVC, Oral/Nasal, Uncuffed	9
Reinforced Tubes	10
Single Use Profile Soft-Seal Cuffed, Murphy Eye, Reinforced Tracheal Tube	10
Single Use Uncuffed, Reinforced Tracheal Tube	10
Bivona™ Aire-Cuf™ Endotracheal Tubes	11
Bivona Uncuffed Endotracheal Tube	11
Cuffed Polar Preformed Tubes	12
Tracheal Tube, Polar™ Preformed, Clear PVC, North Facing, Oral, Profile Soft-Seal (Low Pressure/High Volume) Cuff, Murphy Eye.....	12
Tracheal Tube, Polar Preformed, Ivory PVC, North Facing, Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff.....	12
Tracheal Tube, Polar Preformed, Clear PVC, South Facing, Oral, Profile Soft-Seal (Low Pressure/High Volume) Cuff.....	12
Uncuffed Polar Preformed Tubes	13
Tracheal Tube, Polar Preformed, Clear PVC, North Facing, Oral, Uncuffed, Murphy Eye.....	13
Tracheal Tube, Polar Preformed, Clear PVC, South Facing, Oral, Uncuffed, Murphy Eye	13
Speciality Tubes	14
Endotracheal Tube, Neonatal/Pediatric.....	14
Endotracheal Tube, Neonatal/Pediatric.....	14
Endobronchial Tubes	15
Portex™ Blue Line Endobronchial Tube (Left).....	15
Portex Blue Line Endobronchial Tube (Right)	15
Accessories	16
Portex Nasopharyngeal Airway	16
Portex Intubation Stylet.....	16
Portex Venn Reusable Tracheal Tube Introducer	16
Portex Single Use Bougie Coude.....	16
Portex Tracheal Tube Exchange Guide.....	16
Portex Adult Endotracheal Tube Holder.....	17
Portex Endotracheal Tube Holder (Pediatric)	17



HOME

Cuffed Magill Tubes



Tracheal Tube, Clear PVC, Oral, Murphy Eye, Blue Line, Sacett, Profile Soft-Seal (Low Pressure/High Volume)

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/189/060	6.0	9.1	290	23	10
100/189/065	6.5	9.7	300	23	10
100/189/070	7.0	10.4	310	30	10
100/189/075	7.5	11.3	320	30	10
100/189/080	8.0	11.8	330	30	10
100/189/085	8.5	12.3	330	30	10
100/189/090	9.0	13.0	330	30	10



Tracheal Tube, Clear PVC, Murphy Eye, Oral/Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/199/050	5.0	6.9	250	17	20
100/199/055	5.5	7.6	280	17	20
100/199/060	6.0	8.2	290	23	20
100/199/065	6.5	8.9	300	23	20
100/199/070	7.0	9.6	310	30	20
100/199/075	7.5	10.3	320	30	20
100/199/080	8.0	10.9	330	30	20
100/199/085	8.5	11.6	330	30	20
100/199/090	9.0	12.3	330	30	20
100/199/095	9.5	13.0	330	30	20
100/199/100	10.0	13.7	330	34	20



Cuffed Endotracheal tube with Murphy Eye, Clear PVC, Magill Curve

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/100/030	3.0	4.5	165	10
100/100/035	3.5	5.3	185	10
100/100/040	4.0	6.0	210	10
100/100/045	4.5	6.5	225	10
100/100/050	5.0	6.9	245	10
100/100/055	5.5	7.5	275	10
100/100/060	6.0	8.2	285	10
100/100/065	6.5	8.7	295	10
100/100/070	7.0	9.3	310	10
100/100/075	7.5	10.0	315	10
100/100/080	8.0	10.7	330	10
100/100/085	8.5	11.3	330	10
100/100/090	9.0	12.0	330	10
100/100/010	10.0	13.3	330	10



Cuffed Endotracheal tube with Preloaded Stylet, Clear PVC, Magill Curve, Murphy Eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/102/030	3.0	4.5	165	10
100/102/035	3.5	5.3	185	10
100/102/040	4.0	6.0	210	10
100/102/045	4.5	6.5	225	10
100/102/050	5.0	6.9	245	10
100/102/055	5.5	7.5	275	10
100/102/060	6.0	8.2	285	10
100/102/065	6.5	8.7	295	10
100/102/070	7.0	9.3	310	10
100/102/075	7.5	10.0	315	10
100/102/080	8.0	10.7	330	10
100/102/085	8.5	11.3	330	10
100/102/090	9.0	12.0	330	10
100/102/010	10.0	13.3	330	10




HOME



Tracheal Tube, Siliconized PVC, Oral/Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/166/050	5.0	6.9	250	17	10
100/166/055	5.5	7.6	280	17	10
100/166/060	6.0	8.2	290	23	10
100/166/065	6.5	8.9	300	23	10
100/166/070	7.0	9.6	310	30	10
100/166/075	7.5	10.3	320	30	10
100/166/080	8.0	10.9	330	30	10
100/166/085	8.5	11.6	330	30	10
100/166/090	9.0	12.3	330	30	10
100/166/095	9.5	13.0	330	30	10
100/166/100	10.0	13.7	330	34	10



Tracheal Tube, Ivory PVC, Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/179/050	5.0	7.3	250	17	10
100/179/055	5.5	8.0	280	17	10
100/179/060	6.0	8.8	290	23	10
100/179/065	6.5	9.5	300	23	10
100/179/070	7.0	10.2	310	30	10
100/179/075	7.5	10.9	320	30	10
100/179/080	8.0	11.6	330	30	10
100/179/085	8.5	12.4	330	30	10
100/179/090	9.0	13.1	330	30	10

Uncuffed Magill Tubes



Uncuffed Endotracheal Tube, Murphy Eye, Clear PVC, Magill Curve

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/101/030	3.0	4.5	165	10
100/101/035	3.5	5.3	185	10
100/101/040	4.0	6.0	210	10
100/101/045	4.5	6.5	225	10
100/101/050	5.0	6.9	245	10
100/101/055	5.5	7.5	275	10
100/101/060	6.0	8.2	285	10
100/101/065	6.5	8.7	295	10
100/101/070	7.0	9.3	310	10
100/101/075	7.5	10.0	315	10
100/101/080	8.0	10.7	330	10
100/101/085	8.5	11.3	330	10
100/101/090	9.0	12.0	330	10
100/101/010	10.0	13.3	330	10



Uncuffed Endotracheal Tube, Preloaded Stylet, Clear PVC, Magill Curve, Murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/103/030	3.0	4.5	165	10
100/103/035	3.5	5.3	185	10
100/103/040	4.0	6.0	210	10
100/103/045	4.5	6.5	225	10
100/103/050	5.0	6.9	245	10



Tracheal Tube, Clear PVC, Murphy Eye, Oral/Nasal, Uncuffed				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/127/025	2.5	3.5	150	10
100/127/030	3.0	4.2	170	10
100/127/035	3.5	4.8	190	10
100/127/040	4.0	5.5	210	10
100/127/045	4.5	6.2	230	10
100/127/050	5.0	6.9	250	10
100/127/055	5.5	7.6	280	10
100/127/060	6.0	8.2	290	10
100/127/065	6.5	8.9	300	10
100/127/070	7.0	9.6	310	10
100/127/075	7.5	10.3	320	10
100/127/080	8.0	10.9	330	10
100/127/085	8.5	11.6	330	10
100/127/090	9.0	12.3	330	10
100/127/095	9.5	13.0	330	10
100/127/100	10.0	13.7	330	10



Tracheal Tube, Siliconized PVC, Oral/Nasal, Uncuffed				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/111/020	2.0	3.0	140	10
100/111/025	2.5	3.5	150	10
100/111/030	3.0	4.2	170	10
100/111/035	3.5	4.8	190	10
100/111/040	4.0	5.5	210	10
100/111/045	4.5	6.2	230	10
100/111/050	5.0	6.9	250	10
100/111/055	5.5	7.6	280	10
100/111/060	6.0	8.2	290	10
100/111/065	6.5	8.9	300	10
100/111/070	7.0	9.6	310	10
100/111/075	7.5	10.3	320	10
100/111/080	8.0	10.9	330	10
100/111/085	8.5	11.6	330	10
100/111/090	9.0	12.3	330	10



Tracheal Tube, Siliconized PVC, Murphy Eye, Oral/Nasal, Uncuffed

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/141/025	2.5	3.5	150	10
100/141/030	3.0	4.2	170	10
100/141/035	3.5	4.8	190	10
100/141/040	4.0	5.5	210	10
100/141/045	4.5	6.2	230	10
100/141/050	5.0	6.9	250	10
100/141/055	5.5	7.6	280	10
100/141/060	6.0	8.2	290	10
100/141/065	6.5	8.9	300	10
100/141/070	7.0	9.6	310	10
100/141/075	7.5	10.3	320	10
100/141/080	8.0	10.9	330	10
100/141/085	8.5	11.6	330	10
100/141/090	9.0	12.3	330	10
100/141/095	9.5	13.0	330	10



Tracheal Tube, Siliconized, Oral, Uncuffed

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/126/025	2.5	3.5	113	10
100/126/030	3.0	4.2	123	10
100/126/035	3.5	4.8	133	10
100/126/040	4.0	5.5	143	10
100/126/045	4.5	6.2	153	10
100/126/050	5.0	6.9	163	10
100/126/055	5.5	7.6	173	10
100/126/060	6.0	8.2	193	10



HOME



Tracheal Tube, Ivory PVC, Oral/Nasal, Uncuffed

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/105/025	2.5	3.7	150	10
100/105/030	3.0	4.4	170	10
100/105/035	3.5	5.1	190	10
100/105/040	4.0	5.9	210	10
100/105/045	4.5	6.6	230	10
100/105/050	5.0	7.3	250	10
100/105/055	5.5	8.0	280	10
100/105/060	6.0	8.8	290	10
100/105/065	6.5	9.5	300	10
100/105/070	7.0	10.2	310	10
100/105/075	7.5	10.9	320	10
100/105/080	8.0	11.6	330	10
100/105/085	8.5	12.4	330	10
100/105/090	9.0	13.1	330	10



HOME

Reinforced Tubes



Single Use Profile Soft-Seal Cuffed, Murphy Eye, Reinforced Tracheal Tube

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/110/050	5.0	8.2	259	17	5
100/110/055	5.5	8.6	287	17	5
100/110/060	6.0	9.1	297	23	5
100/110/065	6.5	9.8	307	23	5
100/110/070	7.0	10.4	318	30	5
100/110/075	7.5	11.1	328	30	5
100/110/080	8.0	11.9	338	30	5
100/110/085	8.5	12.4	338	30	5
100/110/090	9.0	13.1	338	30	5
100/110/095	9.5	13.7	338	30	5



Single Use Uncuffed, Reinforced Tracheal Tube

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/113/030	3.0	4.5	175	5
100/113/035	3.5	5.3	195	5
100/113/040	4.0	6.0	215	5
100/113/045	4.5	6.6	245	5
100/113/050	5.0	8.2	259	5



Bivona Cuffed Endotracheal Tubes					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
Bivona TTS Pediatric Endotracheal Tubes					
35W025	2.5	4.0	140	-	1
35W030	3.0	4.7	160	-	1
35W035	3.5	5.4	180	-	1
35W040	4.0	6.0	200	-	1
35W045	4.5	6.6	220	-	1
35W050	5.0	7.3	240	-	1
35W055	5.5	8.0	275	-	1
Bivona Aire-Cuff Endotracheal Tubes					
35W060	6.0	8.7	285	11	1
35W065	6.5	9.4	295	11	1
35W070	7.0	10.0	305	13	1
35W075	7.5	10.4	315	13	1
35W080	8.0	11.0	325	14	1
35W085	8.5	11.8	325	14	1
35W090	9.0	12.3	325	16	1
35W095	9.5	13.3	325	16	1



Bivona Uncuffed Endotracheal Tube				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
15W025	2.5	4.0	140	1
15W030	3.0	4.7	160	1
15W035	3.5	5.4	180	1
15W040	4.0	6.0	200	1
15W045	4.5	6.6	220	1
15W050	5.0	7.3	240	1
15W055	5.5	8.0	275	1
15W060	6.0	8.7	285	1
15W065	6.5	9.4	295	1
15W070	7.0	10.0	305	1
15W075	7.5	10.4	315	1
15W080	8.0	11.0	325	1
15W085	8.5	11.8	325	1
15W090	9.0	12.4	325	1

Cuffed Polar Preformed Tubes



Tracheal Tube, Polar Preformed, Clear PVC, North Facing, Oral, Profile Soft-Seal (Low Pressure/High Volume) Cuff, Murphy Eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/131/060	6.0	8.2	295	23	10
100/131/065	6.5	8.9	305	23	10
100/131/070	7.0	9.6	320	30	10
100/131/075	7.5	10.3	333	30	10
100/131/080	8.0	10.9	345	30	10
100/131/085	8.5	11.6	358	30	10
100/131/090	9.0	12.3	370	30	10



Tracheal Tube, Polar Preformed, Ivory PVC, North Facing, Nasal, Profile Soft-Seal (Low Pressure/High Volume) Cuff

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/133/060	6.0	8.8	410	23	10
100/133/065	6.5	9.5	415	23	10
100/133/070	7.0	10.2	425	30	10
100/133/075	7.5	10.9	435	30	10
100/133/080	8.0	11.6	445	30	10



Tracheal Tube, Polar Preformed, Clear PVC, South Facing, Oral, Profile Soft-Seal (Low Pressure/High Volume) Cuff, Murphy Eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/136/050	5.0	6.9	240	17	10
100/136/055	5.5	7.6	252	17	10
100/136/060	6.0	8.2	265	23	10
100/136/065	6.5	8.9	280	23	10
100/136/070	7.0	9.6	300	30	10
100/136/075	7.5	10.3	315	30	10
100/136/080	8.0	10.9	330	30	10
100/136/085	8.5	11.6	345	30	10
100/136/090	9.0	12.3	360	30	10



HOME

Uncuffed Polar Preformed Tubes



Tracheal Tube, Polar Preformed, Clear PVC, North Facing, Oral, Uncuffed, Murphy Eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/132/030	3.0	4.2	205	10
100/132/035	3.5	4.8	225	10
100/132/040	4.0	5.5	245	10
100/132/045	4.5	6.2	270	10
100/132/050	5.0	6.9	285	10
100/132/055	5.5	7.6	300	10
100/132/060	6.0	8.2	315	10
100/132/065	6.5	8.9	330	10
100/132/070	7.0	9.6	347	10



Tracheal Tube, Polar Preformed, Clear PVC, South Facing, Oral, Uncuffed, Murphy Eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/134/030	3.0	4.2	170	10
100/134/035	3.5	4.8	175	10
100/134/040	4.0	5.5	195	10
100/134/045	4.5	6.2	215	10
100/134/050	5.0	6.9	220	10
100/134/055	5.5	7.6	230	10
100/134/060	6.0	8.2	240	10
100/134/065	6.5	8.9	250	10
100/134/070	7.0	9.6	275	10
100/134/075	7.5	10.3	290	10

Speciality Tubes



Endotracheal Tube, Neonatal/Pediatric				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Tracheal tube				
H4021S	2.5	3.7	155	10
H4022S	3.0	4.4	175	10
H4023S	3.5	4.8	196	10
H4024S	4.0	5.5	216	10
H4025S	4.5	6.1	236	10
H4026S	5.0	6.9	254	10
Tracheal tube with tube holder				
H4003S	3.5	4.8	196	10
H4004S	4.0	5.5	216	10
H4005S	4.5	6.1	236	10
H4006S	5.0	6.9	254	10



Endotracheal Tube, Neonatal/Pediatric				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Tracheal tube				
H4021PS	2.5	3.7	155	10
H4022PS	3.0	4.4	175	10
H4023PS	3.5	4.8	196	10
Tracheal tube with tube holder				
H4003PS	3.5	4.8	196	10

Endobronchial Tubes



Portex Blue Line Endobronchial Tube (Left)

Item Number	Tube Size (Fr)	Tracheal Cuff Resting Diameter (mm)	Bronchial Cuff Resting Diameter (mm)	Bronchoscope Lumen size (mm)	Suction Catheter Size (Fr/mm)	Minimum Equivalent Airway size (mm)	Units per Case
197-28L	28	23	13	3.2	8/2.6	4.6	1
197-32L	32	26	13	3.8	10/3.3	5.4	1
197-35L	35	27	21	4.1	10/3.3	6.0	1
197-37L	37	28	21	4.4	10/3.3	6.3	1
197-39L	39	29	21	4.7	10/3.3	6.6	1
197-41L	41	31	21	5.0	10/4.0	7.0	1

Includes Endobronchial tube with stylet (1), Wye adapter (Carlens style - 1), Suction Catheters (4), Single-axis swivel adapter (1), and 15 mm connectors with auxillary tubes (2).



Portex Blue Line Endobronchial Tube (Right)

Item Number	Tube Size (Fr)	Tracheal Cuff Resting Diameter (mm)	Bronchial Cuff Resting Diameter (mm)	Bronchoscope Lumen size (mm)	Suction Catheter Size (Fr/mm)	Minimum Equivalent Airway size (mm)	Units per Case
198-28R	28	23	13	3.2	8/2.6	4.6	1
198-32R	32	26	13	3.8	10/3.3	5.4	1
198-35R	35	27	21	4.1	10/3.3	6.0	1
198-37R	37	28	21	4.4	10/3.3	6.3	1
198-39R	39	29	21	4.7	10/3.3	6.6	1
198-41R	41	31	21	5.0	10/4.0	7.0	1

Includes Endobronchial tube with stylet (1), Wye adapter (Carlens style - 1), Suction Catheters (4), Single-axis swivel adapter (1), and 15 mm connectors with auxillary tubes (2).



HOME

Accessories




Portex Nasopharyngeal Airway

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
340060	6.0	8.0	146	10
340070	7.0	9.5	168	10
340080	8.0	10.5	177	10
340090	9.0	12.0	187	10



Portex Intubation Stylet

Item Number	Recommended Tube I.D. (mm)	Size (Fr/mm)	Length (in/mm)	Units per Case
103004	2.5 - 4.5	6/2.0	11.75/298	10
103005	5.0 - 7.0	8/2.6	13.75/349	10
103006	7.5 - 10.0	10.5/3.5	15.75/400	10




Portex Venn Reusable Tracheal Tube Introducer

Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (in/mm)	Units per Case
14-504-17	6.0 – 11.0	5.0	23.62/600	1



Portex Single Use Bougie Coude

Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (in/mm)	Units per Case
100/123/515	6.0 – 11.0	5.0	27.55/700	10



Portex Tracheal Tube Exchange Guide

Item Number	Recommended Tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable				
14-504-76	4.0 – 6.0	3.3	700	1
14-504-68	6.0 – 11.0	5.0	700	1
Single-use				
100/125/005	2.5 – 3.5	1.7	500	10
100/125/010	4.0 – 6.0	3.3	700	10
100/123/015	6.0 – 11.0	5.0	700	10



Portex Adult Endotracheal Tube Holder		
Item Number	Item Description	Units per Case
K1000	Holder for 4.5 – 6.5 mm	30
K2000	Holder for 7.0 – 8.5 mm	30
K3000	Holder for 9.0 – 10.0 mm	30



Portex Endotracheal Tube Holder (Pediatric)		
Item Number	Item Description	Units per Case
100/603/025	2.5 mm Tracheal Tube Holder	10
100/603/030	3.0 mm Tracheal Tube Holder	10
100/603/035	3.5 mm Tracheal Tube Holder	10
100/603/040	4.0 mm Tracheal Tube Holder	10
100/603/045	4.5 mm Tracheal Tube Holder	10

Index

100/100/010.....4	100/102/065.....4	100/111/025.....7	100/127/050.....7
100/100/030.....4	100/102/070.....4	100/111/030.....7	100/127/055.....7
100/100/035.....4	100/102/075.....4	100/111/035.....7	100/127/060.....7
100/100/040.....4	100/102/080.....4	100/111/040.....7	100/127/065.....7
100/100/045.....4	100/102/085.....4	100/111/045.....7	100/127/070.....7
100/100/050.....4	100/102/090.....4	100/111/050.....7	100/127/075.....7
100/100/055.....4	100/103/030.....6	100/111/055.....7	100/127/080.....7
100/100/060.....4	100/103/035.....6	100/111/060.....7	100/127/085.....7
100/100/065.....4	100/103/040.....6	100/111/065.....7	100/127/090.....7
100/100/070.....4	100/103/045.....6	100/111/070.....7	100/127/095.....7
100/100/075.....4	100/103/050.....6	100/111/075.....7	100/127/100.....7
100/100/080.....4	100/105/025.....9	100/111/080.....7	100/131/060.....12
100/100/085.....4	100/105/030.....9	100/111/085.....7	100/131/065.....12
100/100/090.....4	100/105/035.....9	100/111/090.....7	100/131/070.....12
100/101/010.....6	100/105/040.....9	100/113/030.....10	100/131/075.....12
100/101/030.....6	100/105/045.....9	100/113/035.....10	100/131/080.....12
100/101/035.....6	100/105/050.....9	100/113/040.....10	100/131/085.....12
100/101/040.....6	100/105/055.....9	100/113/045.....10	100/131/090.....12
100/101/045.....6	100/105/060.....9	100/113/050.....10	100/132/030.....13
100/101/050.....6	100/105/065.....9	100/123/015.....16	100/132/035.....13
100/101/055.....6	100/105/070.....9	100/123/515.....16	100/132/040.....13
100/101/060.....6	100/105/075.....9	100/125/005.....16	100/132/045.....13
100/101/065.....6	100/105/080.....9	100/125/010.....16	100/132/050.....13
100/101/070.....6	100/105/085.....9	100/126/025.....8	100/132/055.....13
100/101/075.....6	100/105/090.....9	100/126/030.....8	100/132/060.....13
100/101/080.....6	100/110/050.....10	100/126/035.....8	100/132/065.....13
100/101/085.....6	100/110/055.....10	100/126/040.....8	100/132/070.....13
100/101/090.....6	100/110/060.....10	100/126/045.....8	100/133/060.....12
100/102/010.....4	100/110/065.....10	100/126/050.....8	100/133/065.....12
100/102/030.....4	100/110/070.....10	100/126/055.....8	100/133/070.....12
100/102/035.....4	100/110/075.....10	100/126/060.....8	100/133/075.....12
100/102/040.....4	100/110/080.....10	100/127/025.....7	100/133/080.....12
100/102/045.....4	100/110/085.....10	100/127/030.....7	100/134/030.....13
100/102/050.....4	100/110/090.....10	100/127/035.....7	100/134/035.....13
100/102/055.....4	100/110/095.....10	100/127/040.....7	100/134/040.....13
100/102/060.....4	100/111/020.....7	100/127/045.....7	100/134/045.....13



HOME

100/134/050.....13	100/166/080.....5	100/603/045.....17	340090.....16
100/134/055.....13	100/166/085.....5	103004.....16	35W025.....11
100/134/060.....13	100/166/090.....5	103005.....16	35W030.....11
100/134/065.....13	100/166/095.....5	103006.....16	35W035.....11
100/134/070.....13	100/166/100.....5	14-504-17.....16	35W040.....11
100/134/075.....13	100/179/050.....5	14-504-76.....16	35W045.....11
100/136/050.....12	100/179/055.....5	14-504-68.....16	35W050.....11
100/136/055.....12	100/179/060.....5	15W025.....11	35W055.....11
100/136/060.....12	100/179/065.....5	15W030.....11	35W060.....11
100/136/065.....12	100/179/070.....5	15W035.....11	35W065.....11
100/136/070.....12	100/179/075.....5	15W040.....11	35W070.....11
100/136/075.....12	100/179/080.....5	15W045.....11	35W075.....11
100/136/080.....12	100/179/085.....5	15W050.....11	35W080.....11
100/136/085.....12	100/179/090.....5	15W055.....11	35W085.....11
100/136/090.....12	100/189/060.....3	15W060.....11	35W090.....11
100/141/025.....8	100/189/065.....3	15W065.....11	35W095.....11
100/141/030.....8	100/189/070.....3	15W070.....11	H4003PS.....14
100/141/035.....8	100/189/075.....3	15W075.....11	H4003S.....14
100/141/040.....8	100/189/080.....3	15W080.....11	H4004S.....14
100/141/045.....8	100/189/085.....3	15W085.....11	H4005S.....14
100/141/050.....8	100/189/090.....3	15W090.....11	H4006S.....14
100/141/055.....8	100/199/050.....3	197-28L.....15	H4021PS.....14
100/141/060.....8	100/199/055.....3	197-32L.....15	H4021S.....14
100/141/065.....8	100/199/060.....3	197-35L.....15	H4022PS.....14
100/141/070.....8	100/199/065.....3	197-37L.....15	H4022S.....14
100/141/075.....8	100/199/070.....3	197-39L.....15	H4023PS.....14
100/141/080.....8	100/199/075.....3	197-41L.....15	H4023S.....14
100/141/085.....8	100/199/080.....3	198-28R.....15	H4024S.....14
100/141/090.....8	100/199/085.....3	198-32R.....15	H4025S.....14
100/141/095.....8	100/199/090.....3	198-35R.....15	H4026S.....14
100/166/050.....5	100/199/095.....3	198-37R.....15	K1000.....17
100/166/055.....5	100/199/100.....3	198-39R.....15	K2000.....17
100/166/060.....5	100/603/025.....17	198-41R.....15	K3000.....17
100/166/065.....5	100/603/030.....17	340060.....16	
100/166/070.....5	100/603/035.....17	340070.....16	
100/166/075.....5	100/603/040.....17	340080.....16	

