



NONIN

Nonin Pulse Oximeter Sensor Compatibility Guide

Use the chart below to find the right sensor for your needs. Visit nonin.com for more information and video tutorials.




	MODEL NUMBER	DESCRIPTION	CABLE LENGTH	APPLICATION	PACKAGING	WristOx ₂ ® Model 3150	CO-Pilot™*	PalmSAT® 2500/2500A	8500 Series	Model 7500**	Avant® 9700
--	--------------	-------------	--------------	-------------	-----------	-----------------------------------	------------	---------------------	-------------	--------------	-------------

Reusable Soft Sensors **8000S Series**

	8000SS 8000SS-3	Small Soft Sensor 7.5 – 12.5 mm digit 0.3 – 0.5 inches digit	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
	8000SS-WO2	Small Soft Sensor – WristOx₂ 7.5 – 12.5 mm digit 0.3 – 0.5 inches digit	0.3m	Index, middle, or ring finger	Single	✓					
	8000SM 8000SM-3	Medium Soft Sensor 10 – 19 mm digit 0.4 – 0.75 inches digit	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
	8000SM-WO2	Medium Soft Sensor – WristOx₂ 10 – 19 mm digit 0.4 – 0.75 inches digit	0.3m	Index, middle, or ring finger	Single	✓					
	8000SL 8000SL-3	Large Soft Sensor 12.5 – 25.5 mm digit 0.5 – 1.0 inches digit	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
	8000SL-WO2	Large Soft Sensor – WristOx₂ 12.5 – 25.5 mm digit 0.5 – 1.0 inches digit	0.3m	Index, middle, or ring finger	Single	✓					

Reusable Finger Clip, Ear Clip and Reflective Sensors **8000 Series**

	8330AA	Finger Clip Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle, or ring finger	Single		✓				
	8000AA 8000AA-2 8000AA-3	Finger Clip Sensor – Adult > 30 kg / > 66 lbs	1m 2m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
	8000AA-WO2	Finger Clip Sensor – WristOx₂ > 30 kg / > 66 lbs	0.3m	Index, middle, or ring finger	Single	✓					
	8000AP 8000AP-3	Finger Clip Sensor – Pediatric 8 – 30 kg / 18 – 66 lbs	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
	8000Q2 8000Q2-WO2	Ear Clip Sensor > 40 kg / > 88 lbs	1m	Earlobe	Single			✓	✓	✓	✓
	8000R 8000H	Reflectance Sensor > 30 kg / > 66 lbs Sensor Holder 10 caps /20 stickers	1m N/A	Forehead	Single			✓	✓	✓	✓

Reusable Flex Sensors with FlexiWrap **8000 Series**

	8000J 8000J-3 8000JFW	Reusable FlexSensor – Adult > 22 kg / > 44 lbs incl. 25 Adult FlexiWraps	1m 3m N/A	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
	8000J-WO2 8000JFW	Reusable FlexSensor – WristOx₂ > 22 kg / > 44 lbs incl. 25 Adult FlexiWraps	0.3m N/A	Index, middle, or ring finger	Single 24 per pack	✓					
	8008J 8008JFW	Reusable FlexSensor – Infant 2 – 22 kg / 4.4 – 44 lbs incl. 25 Infant FlexiWraps	1m N/A	Toe	Single 24 per pack			✓	✓	✓	✓
	8001J 8001JFW	Reusable FlexSensor – Neonatal > 2 kg / 4.4 lbs incl. 25 Neonatal FlexiWraps	1m N/A	Side of foot	Single 24 per pack			✓	✓	✓	✓

Single-Patient-Use Cloth Sensors **6000C Series**

	6000CA	Cloth Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
	6000CA-WO2	Cloth Sensor – Adult > 30 kg / > 66 lbs	0.3m	Index, middle, or ring finger	Single 15 per pack	✓					
	6000CP	Cloth Sensor – Pediatric 10 – 40 kg / 22 – 88 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
	6000CI	Cloth Sensor – Infant 2 – 20 kg / 4.4 – 44 lbs	1m	Big toe	Single 24 per pack			✓	✓	✓	✓
	6000CN	Cloth Sensor – Neonatal/Adult Neonatal: < 2 kg / < 4.4 lbs Adult: > 30 kg / > 66 lbs	1m	Foot	Single 24 per pack			✓	✓	✓	✓

Single-Patient-Use Flexi-Form® III Sensors **7000 Series**

	7000A	Flexi-Form Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
	7000P	Flexi-Form Sensor – Pediatric 10 – 40 kg / 22 – 88 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
	7000I	Flexi-Form Sensor – Infant 2 – 20 kg / 4.4 – 44 lbs	1m	Big toe	Single 24 per pack			✓	✓	✓	✓
	7000N	Flexi-Form Sensor – Neonatal/Adult Neonatal: < 2 kg / < 4.4 lbs Adult: > 30 kg / > 66 lbs	1m	Foot	Single 24 per pack			✓	✓	✓	✓

Accessories

	UNI-EXT-1 UNI-EXT-3 UNI-EXT-6 UNI-EXT-9	Extension Cable	1m 3m 6m 9m		Single			✓	✓	✓	✓
	UNI-RA-3 UNI-RA-0	Right Angle Connector	3m 15cm		Single			✓	✓	✓	✓
	3150I	Adapter Cable to 9-pin (D-Connector). Connects PureLight® sensors to WristOx ₂ Model 3150	12cm		Single	✓					



*Also available is the SenSmart Model 8204CA Regional Oximetry Sensor for use with the Nonin CO-Pilot.
**Also available is the Model 7500FO for use in MRI environments utilizing the 8000FI and 8000FC cable and sensor.
© 2022 Nonin Medical, Inc. All trademarks are the property of Nonin Medical, Inc. unless otherwise noted.
M-19-028-03