



FLUSSO™

FLUSSO™ BY PASS ADAPTER

To Maintain a Closed Circuit



Integrated Swing Valve Technology™



www.keepthepeep.com



FLUSSO™ BY PASS ADAPTER

Flusso™ facilitates the safe disconnection of the patient from the mechanical ventilator during a circuit change, HME change, or disconnection of patient for transport.

Flusso™ By Pass allows for the maintenance of a patient's positive end expiratory pressure (PEEP) and reduces potential exhaled aerosols during ventilator circuit disconnect.

ORDERING INFORMATION

| Item # | Description | Qty |
|----------|-------------------------|---------|
| 420-1978 | Flusso™ By Pass Adapter | 25/case |

- Proprietary valve with Swing Valve Technology™
- Clear housing with coloured valve permits visualization of function of Swing Valve Technology™
- Integrated swivel at patient connection to reduce torque on endotracheal tube
- Tethered dust cap keeps ports clean and protects swing valve
- Minimal dead space (9cc)
- Standard ISO 15/22mm connections

Maintain PEEP • Reduce Exhaled Aerosols • Lung Protective Ventilation



McArthur Medical Sales Inc
1846 5th Concession Rd W
Rockton, ON L0R 1X0

1-800-996-6674
mmsi@mcarthurmedical.com
www.mcarthurmedical.com

FLUSSO™
www.keepthepeep.com